

FISSATIVI - *FIXATIVES*

Formaldeide 40% *Formaldehyde 40%*

CE IVD

Cod.	Descrizione Description	Q.ta Q.ty
1625	stabilizzata bottiglia/ <i>stabilized bottle</i>	1 lt.
1626	stabilizzata tanica/ <i>stabilized tank</i>	5 lt.
1620	stabilizzata tanica/ <i>stabilized tank</i>	25 lt.
1628	tamponata tanica/ <i>buffered tank</i>	5 lt.
1629	tamponata tanica/ <i>buffered tank</i>	10 lt.



Formaldeide soluzione 8% *Formaldehyde 8%*

CE IVD

Cod.	Descrizione Description	Q.ta Q.ty
1612	bottiglia/ <i>bottle</i>	1 lt.
1613	tanica/ <i>tank</i>	5 lt.
1618	tanica/ <i>tank</i>	10 lt.
1619	tanica/ <i>tank</i>	20 lt.

Formaldeide soluzione 4% pH 6,9-7,1 *Formaldehyde 4% pH 6,9-7,1*

CE IVD

Cod.	Descrizione Description	Q.ta Q.ty
1609	bottiglia/ <i>bottle</i>	1 lt.
1614	tanica/ <i>tank</i>	5 lt.
1617	tanica/ <i>tank</i>	10 lt.
1616	tanica/ <i>tank</i>	20 lt.

Formaldeide soluzione 4% pH 7,2-7,4 *Formaldehyde 4% pH 7,2-7,4*

CE IVD

Cod.	Descrizione Description	Q.ta Q.ty
1636	tanica/ <i>tank</i>	5 lt.
1637	tanica/ <i>tank</i>	10 lt.
1638	tanica/ <i>tank</i>	20 lt.

Cod.	Descrizione Description	Q.ta Q.ty
3626	rubinetto erogatore per tanica/ <i>tap for tank</i>	1 pcs

