

\ Treatment \ Post-operative Redon systems \ Low vacuum system  
**500ml**



**low vacuum system for routine use in operating theatre**

product features:

- starting with a vacuum of 120 mbar
- content: 500 ml
- connecting tube, length: 125 cm
- universal stepped connector for safe connection of ch. 6 to ch. 18 redon drains
- set with drain and guiding needle
- length of the drain: 800 mm
- perforation of the drain: 150 mm
- individually packed and sterilized
- special designs available on request
- UMDNS-Code: 11-307
- Klasse/Regel: IIa/7

**Art.-Nr.: 21599, PRI-LOW-VAC Redon set, 500 ml with drain (800/150) ch. 14 and guiding needle**

Packing unit (VPE): 35, Shipping unit (VSE): 35



**low vacuum system for routine use in operating theatre**

product features:

- starting with a vacuum of 120 mbar
- content: 500 ml
- connecting tube, length: 125 cm
- universal stepped connector for safe connection of ch. 6 to ch. 18 redon drains
- set with drain and guiding needle
- length of the drain: 800 mm
- perforation of the drain: 150 mm
- individually packed and sterilized
- special designs available on request
- UMDNS-Code: 11-307
- Klasse/Regel: IIa/7

**Art.-Nr.: 21598, PRI-LOW-VAC Redon set, 500 ml with drain (800/150) ch. 16 and guiding needle**

Packing unit (VPE): 35, Shipping unit (VSE): 35