

12 tests (Réf. 44412)

Description	Quantité
SUSPENSION FUNGI : Flacon de 4 mL de milieu semi-gélosé tamponné pour la mise en suspension des colonies et la standardisation; contenant du Bacto Agar, de la colimycine et de la vancomycine.	14
MES FUNGI : Flacon de 5 mL de milieu liquide d'étude de la sensibilité (milieu RPMI modifié) contenant de la résazurine comme indicateur de croissance et du glucose.	12
TC FUNGI : Flacon de 4 mL de sulfate de baryum pour le contrôle de standardisation de l'inoculum.	1
FUNGIFAST AFG : Galerie de 20 puits, pour un test, conditionnée individuellement en sachet aluminium.	12

12 tests (Ref. 44412)

Description	Quantity
SUSPENSION FUNGI : 4 mL-vial of buffered semi-agar medium containing Bacto Agar, colimycin and vancomycin, for colony suspension and standardization.	14
MES FUNGI : 5 mL-vial of liquid medium (modified RPMI medium) for susceptibility testing, containing resazurin (growth indicator) and glucose.	12
TC FUNGI : 4 mL vial of barium sulphate solution for the control of inoculum standardization.	1
FUNGIFAST AFG : 20 well-tray, allowing the testing of one sample, individually packed in an aluminium sachet	12

Parc d'activités du Plateau
19, allée d'Athènes
83870 SIGNES
FRANCE
Tél. : 33 (0)4 94 88 55 00
Fax : 33 (0)4 94 88 55 22

