

12. Mikrokategoriai

Headway 156cm
Microcatheter

Headway®
Now with Duo! Technology
microcatheter

Hybrid braid/coil construction to improve stent delivery experience



Headway 21 156cm
Microcatheter

Headway 27 156cm
Microcatheter

LVIS
Intraluminal Support Device

ERIC
Endovascular Retractor with Irrigated Cages

FRED Jr.
Flow Re-Direction Endovascular Device

Sofia
Distal Access Catheter
Soft torqueable catheter Optimized For Intracranial Access

FRED
Flow Re-Direction Endovascular Device

Sofia
Distal Access Catheter
Soft torqueable catheter Optimized For Intracranial Access

Stretch Resistant Shaft

- For stability and prevention of catheter body twisting

Hybrid Braid & Coil Design

- Coil reinforcement provides lumen integrity, bending flexibility and excellent shape retention
- Proximal variable braid reinforcement provides support and torque control
- Hybrid design provides 1:1 push/pull control for trackability

Hibridinė pynės ir spiralės konstrukcija

Soft Distal Shaft

- Smooth and atraumatic tracking

PTFE Liner

- Lubricious, durable inner liner for smooth stent delivery

Dėl sienelėje siaurėjančio PTFE arba lygiavertės* medžiagos sluoksnio, sustiprinto spirale, galima 1:1 "push/pull" kateterio kontrolė; Sienelėje esanti spirale neleidžia kateteriui ovalizuotis ir užtikrina maksimalų lankstumą;

Lubricious Hydrophilic Coating

- Reduced friction during navigation in tortuous anatomy and during lengthy procedures

Paviršiaus hidrofilinės dangos ilgis – 100 cm.

SPECIFICATIONS

Headway® 21 & 27 (156cm) Microcatheter

Microcatheter with Hydrophilic Coating

Packed 1 per box; includes shaping mandrel and introducer sheath

Product Name	Description	Product Code	ID (inches)	Tip Shape	Usable Length (cm)	OD Prox./Distal (French)	Tip Markers
Headway 21 156cm	Straight	MC212156S	0.021	STR	156	2.5 / 2.0	2
Headway 27 156cm	Straight	MC272156S	0.027	STR	156	3.1 / 2.6	2

Mikrokategorio vidinis skersmuo – 0,0165"; 0,017"; 0,021"; 0,027".


Proksimalinis/distalinis skersmuo OD: 2,1Fr/1,3 Fr; 2,1 Fr/1,6 Fr; 2,4 Fr/1,7 Fr; 2,5 Fr/2,0 Fr; 3,1 Fr/2,6 Fr.

Bendras mikrokategorio ilgis (siūlyti ne mažiau kaip 2 ilgių): 156 ±3 cm ir 167 ±3 cm.



Paviršiaus hidrofilinės dangos ilgis – 100 cm

turi 2 mikrokategorio galiuko žymėjimus.



Hybrid Braid/Coil Design

Features

- Proximal column strength for stability and torqueability
- 156cm working length with distal flexible transitions
- Resistance to kinking and ovalization
- Resistance to stretching
- Resistance to twisting

Benefits

- ✓ Facilitate access to the lesion
- ✓ Smooth device trackability
- ✓ Controlled device delivery

INDICATIONS FOR USE

The Headway Microcatheter is intended for general intravascular use, including the peripheral, coronary and neuro vasculature for the infusion of diagnostic agents, such as contrast media, and therapeutic agents, such as occlusion coils.

MICROVENTION, Headway, LVIS, FRED, Sofia and ERIC are registered trademarks of MicroVention, Inc.
 • Scientific and clinical data related to this document are on file at MicroVention, Inc. Refer to Instructions for Use, contraindications and warnings for additional information. ©2016 MicroVention, Inc. MM457(i) 3/16



MicroVention, Inc.
Worldwide Headquarters
 1311 Valencia Avenue
 Tustin, CA 92780 USA
 MicroVention UK Limited
 MicroVention Europe, S.A.R.L.
 MicroVention Deutschland GmbH
 Web

PH +1.714.247.8000

PH +44 (0) 191 258 6777
 PH +33 (1) 39 21 77 46
 PH +49 211 210 798-0
 microvention.com