

## 8. Intermediate kateteriai



# Distal Access Catheter

**Enhanced control and stability for superior microcatheter performance and reliable device delivery**

**SOFIA™**

Distal Access Catheter

Soft Torqueable Catheter Optimized For Intracranial Access

1 per box and includes shaping mandrel and introducer sheath

Product Code	Catheter Size (French)	Proximal OD (inch/mm)	Distal OD (inch/mm)	ID (inch)	Working Length (cm)	Distal Length (cm)	Proximal Length (cm)	Distal Tip Shape
DA5115ST	5F	0.068 / 1.7	0.068 / 1.7	0.055	115	17	98	STRAIGHT
DA5125ST	5F	0.068 / 1.7	0.068 / 1.7	0.055	125	17	108	STRAIGHT

Distaliai matmenys – OD ne mažiau 0,068 / 1,7  
Proksimaliniai matmenys – OD ne mažiau 0,068 / 1,7;

Vidinis  
diametras  
ID 0,055"

Proksimalaus galo ilgis: 96 cm, 98 cm  
kai kateterio darbinis ilgis 115 cm; 108 cm, 106 cm  
kai kateterio darbinis ilgis 125 cm; 112 cm  
kai kateterio darbinis ilgis 131 cm.  
Darbinis kateterio ilgis: 115 cm; 125 cm; 131cm;

## Features

Exceptionally Soft Distal Shaft

Low Profile

Steam Shapeable Tip and  
Torqueable Shaft

Enhanced Kink Resistance

## Benefits

Easier navigation in tortuous vessels

Compatible with 6F .070 inch ID guide catheter systems

Steerable around bifurcations

Easy to handle with 1:1 push/pull control

Spiralės konstrukcija turi užtikrinti  
stabilumą 1:1 stumiant ir traukiant;



MicroVention, Inc.

Worldwide Headquarters

1311 Valencia Avenue

Tustin, CA 92780, USA

MicroVention UK Limited

MicroVention Europe, S.A.R.L.

MicroVention Deutschland GmbH

Web

PH +1.714.247.8000

PH +44 (0) 191 258 6777

PH +33 (1) 39 21 77 46

PH +49 211 210 798-0

microvention.com



MICROVENTION is a registered trademark of MicroVention, Inc. SOFIA is a trademark of MicroVention, Inc.

• Scientific and clinical data related to this document are on file at MicroVention, Inc. Refer to Instructions for Use for additional information. ©2014 MicroVention, Inc. MM323(I) Rev.A 9/14