

Eliminate™ components

- x2 Locking syringes 30 ml polycarbonate
- x1 Extension line with 2-way stopcock
- x1 Flushing tool
- x2 Filter baskets 70 µm

Ordering information

Guide catheter compatibility

Eliminate™ with stylet
6Fr (I.D. ≥ 0.70 inch / 1.78 mm)
7Fr (I.D. ≥ 0.80 inch / 2.03 mm)
EG1602
EG1652

Product specifications

Tip Length	6 mm
Coating	Hydrophilic
Rapid exchange segment	23 cm
Radiopaque marker band	1 mm located at 4 mm from the tip
Positioning marks	10 cm = single white mark located at 90 cm from distal tip
Wire compatibility	0.014" (0.36 mm)
Catheter length	140 cm

TERUMO CORPORATION
44-1, 2-chome
Hongo
Tokyo 113-8501
Japan
Tel. +81-3-5374-8111
Fax +81-3-5374-8196

European Headquarters and Manufacturing
TERUMO EUROPE N.V.
3001 Leuven
Belgium
Tel. +32-16-40122-49
Fax +32-16-40122-49

European Sales Offices
BENELUX SALES DIVISION
3001 Leuven
Belgium
Tel. +32-16-40122-49
Fax +32-16-40122-49

LABORATOIRES TERUMO FRANCE S.A.
Place Armand - 834, Urmess
92100 Nanterre
France
Tel. +33-1-4700-1100
Fax +33-1-4700-1100

TERUMO DEUTSCHLAND GMBH
Zugspitzenstrasse 10
D-42699 Solingen
Germany
Tel. +49-212-374-1100
Fax +49-212-374-1100

TERUMO EUROPE SPAIN SL
Avda. Juan Carlos I
28002 Madrid
Spain
Tel. +34-91-555-1100
Fax +34-91-555-1100

TERUMO ITALIA SRL
Via Vercelli, 10
I-20134 Milano
Italy
Tel. +39-02-5741-1100
Fax +39-02-5741-1100

TERUMO DEUTSCHLAND GMBH
Hauptstrasse 87
D-42699 Solingen
Germany
Tel. +49-212-374-1100
Fax +49-212-374-1100

TERUMO SWEDEN AB
Nya Värmdö, Björnså 60
SE-425 71 Mölndal
Sweden
Tel. +46-31-748-5800
Fax +46-31-748-5800

TERUMO UK LTD
The Priory
Brynmor, Surrey
Surrey, Surrey
UK
Tel. +44-1753-476-200
Fax +44-1753-476-200

TERUMO EUROPE N.V.
Mitsubishi Representative Office
121108 Moscow
Russia
Tel. +7-495-344-3300
Fax +7-495-344-3303

TERUMO
Office in Russia
Mitsubishi Representative Office
121108 Moscow
Russia
Tel. +7-495-344-3300
Fax +7-495-344-3303

NY Regional Representative of Terumo Corporation - Published by Terumo Europe N.V.
Address: ny@terumo.com
Change: please call our website.

ELIMINATE™

Aspiration catheter

2" 15 7 4 3 8 0

► The deliverable solution for thrombus removal ◀




The deliverable solution for thrombus removal

Eliminate™ is an innovative aspiration catheter designed to offer an **optimal balance** between **crossing** performance, **kink resistance** and **thrombo-aspiration capability**

>>> Deliverability

Eliminate's innovative technology provides improved deliverability and kink resistance

- **Pre-loaded stylet** provides improved pushability and ensures optimal crossability in difficult to cross scenarios
- **Fully braided shaft construction** for improved kink resistance of the catheter
- **23 cm long rapid exchange segment** provides improved trackability
- **Hydrophilic coating** facilitates easy navigation through tortuous anatomies



...**ELIMINATE™**
ASPIRATION CATHETER

6 Fr guide catheter compatibility




7 Fr guide catheter compatibility



>>> Aspiration Capability

Advanced technology provides enhanced aspiration

- **Large extraction lumen** and the braided shaft construction ensures constant aspiration performance throughout the procedure
- **Rounded and short tip design** preserves crossability, while providing safe and effective aspiration
- A choice of **2 sizes** (6 Fr or 7 Fr guide catheter compatibility) to meet different clinical scenario



>>> Ease of Use

Single-handed operator

- **Short and thin distal radiopaque marker** provides an excellent distal visualization and increases distal flexibility
- **Complete kit contains:**
 - 30 ml locking syringes
 - Extension line (25 cm) with 2 way stopcock
 - Filters basket

