

# CKG16, Calibration kit


## DESCRIPTION

**CKG16** for glass tubes with diameter 16 mm, set of 0.5; 1.0; 2.0; 3.0; 4.0 McFarland Turbidity Standards (latex particles).

## SPECIFICATIONS

Diameter

16 mm



## CAT. NUMBER

BS-050102-BK

CKG16 - ø16mm - 0.5; 1.0; 2.0; 3.0; 4.0 McF

## ACCESSORIES



### DEN-1

BS-050102-AAF

Densitometer (suspension  
turbidity detector)

Densitometer is designed for measurement of cell suspension's turbidity in the range of 0.0–6.0 McFarland units ( $0 - 180 \times 10^7$  cells/ml).

Densitometer provides the opportunity to measure solution ...

[read more](#)



### DEN-1B

BS-050104-AAF

Densitometer (suspension  
turbidity detector)

Densitometer is designed for measurement of cell suspension's turbidity in the range of 0.0–6.0 McFarland units ( $0 - 180 \times 10^7$  cells/ml).

Densitometers provide the opportunity to measure solution turbidity in ...

[read more](#)



### DEN-600

BS-050109-AAA

Photometer

**DEN-600** is a compact, portable, rechargeable battery powered photometer. DEN-600 comprises of 600 nm wavelength optical system, which enables to apply - 1) OD600 method that estimated total number of ...

[read more](#)



### DEN-600m

IVD certified

Photometer IVD

**DEN-600m** (IVD version) is a compact, portable, rechargeable battery powered photometer. DEN-600m comprises of 600 nm wavelength optical system, which enables to apply - 1) OD600 method that estimated total ...

[read more](#)