

# tg<sup>®</sup> fix

## Tubular Net Bandage



### Fields of application

- for retention of all types of wound/ointment dressings
- for application of occlusion films
- for retention of splints

### Properties

- wide-mesh knitted fabric
- bi-elastic
- does virtually not crease when adjusted
- less tearing after being cut to size

### Product composition

71% polyamide (crinkled), 29% non-brightened elastodiene (latex)

### To note

- washable
- can be sterilised with water vapour (121°C) and radiation

## tg<sup>®</sup> fix

stretched length

Size	REF	Shipping units
<b>4 m rolled in a folding box</b>		
A for finger	24240	1/100
B for small extremities (several fingers together/hand/foot)	24241	1/100
C for large extremities (small head/arm/leg)	24242	1/1/98
D for large head, small trunk	24243	1/58
E for large trunk, hip, armpit	24244	1/1/28
<b>25 m folded in a cardboard dispenser</b>		
A for finger	24250	1/66
B for small extremities (several fingers together/hand/foot)	24251	1/20
C for large extremities (small head/arm/leg)	24252	1/1/16
D for large head, small trunk	24253	1/16
E for large trunk, hip, armpit	24254	1/12

Not available in all countries