



Huber needles range with and without connecting line

Huber needles

- Complete range of sizes and gauges
- For the puncture of implantable ports

Description	Gauge	Needle Ø (mm)	Usable length (mm)	Product codes
Curved	23	0.6	20	522506
			25	522507
	22	0.7	35	523807
			20	522509
	20	0.9	35	523809
			20	522511
Straight	23	0.6	35	523811
			20	522506
	22	0.7	30	512507
			38	513807
	20	0.9	25	512509
			25	512511
	19	1.1	28	513811



polyperf™ without Y Site - Huber needle

- Curved Huber needle with connecting line (25 cm)
- No interaction with anti-mitotic drugs

Gauge	Needle Ø (mm)	Usable length (mm)	Product codes
22	0.7	15	581507
		17	581707
		20	582007
		25	582507
		30	583007
		35	583507
20	0.9	15	581509
		17	581709
		20	582009
		25	582509
		30	583009
		35	583509
19	1.1	15	581511
		17	581711
		20	582011
		25	582511
		20	583011
		35	583511

Winglets:

- good grip
- good stability



Winglets colour coded:

- easy identification of the gauge

Bold references are the most commonly used. Other references are only available upon request with a minimum lead time of 12 weeks.

polyperf™ with Y site - Huber needle

Gauge	Needle Ø (mm)	Usable length (mm)	Product codes
22	0.7	20	592007
		25	592507
20	0.9	20	592009
		25	592509
19	1.1	20	592011
		25	592511



Packed in box of 10 units

Bold references are the most commonly used. Other references are only available upon request with a minimum lead time of 12 weeks.

This device is not made with dry or natural rubber latex.

Huber needles and polyperf™ are sterilized by ethylene oxide. polyperf™ is a registered trademark of PEROUSE MEDICAL.