


Universal chair FK-02 is intended for injection procedures. The frame of the medical procedure chair is made of metal structures, painted with powder polymer paint and has four height-adjustable legs. (8.2) Two ball-bearing adjustable armrests make it easier to carry out the injection procedure, both on the left and right side of the patient. (8.7) The backrest of the chair is adjustable with the help of a gas spring. (8.1) All parts of the chair are made of health-friendly, non-allergenic materials, easy to clean and disinfect, resistant to fire and mold (8.4). The legs of the chair are covered with transparent plastic that is resistant to friction and disinfectants. (8.3).

Parameters:	Model FK02
Back support adjustment angle (8.8)	0° - 95°
Overall length (8.6)	170 cm
Seat width (8.6)	50cm
Total width of the chair (8.6)	71 cm
Trendelenburg position (8.9)	from 0° to -3°
Handrail adjustment using a ball bearing (8.7)	from -15° to+ 65°
Permissible load (8.10)	150 kg



R. 6  
•1

••

•••


WST

w R



---

WSTECHS.C. R. Lorens, W. Pielech, J. Mazur, 38- 500 Sanok ul. Okulickiego 43



**{8.5}** Chair table 25x35 cm; load 12 kg