



STATEMENT

Ravimed Sp. z o.o. confirms that:

- Each set is packed individually into a single, transparent bag. Then individually packed 3 sets into one collective bag.
- The tube of erythrocyte storage bag is 600mm, individually numbered every 40mm.
- Ravimed blood bags are manufactured in accordance to ISO 3826 standards. So the bags have rounded corners to prevent friction during separation and transfusion of blood components and have 2 slits on both sides of all blood bags in the system.
- Ravimed blood bags has 3 years of expiry date and the quality certificate is made for all the life time of the product. At the same time Ravimed recommends to store at a room temperature. When you store at a room temperature - the anticoagulant is colorless all 3 years.
- Ravimed needle meets the requirements mentioned below:
 - Needle is made of stainless steel thin walled siliconized 16 G, ensuring sufficient blood flow;
 - The needle cover reliably ensures the hermeticity of the needle without the possibility of self-sealing;
 - Convenient picking up of the needle holder ensures high quality venipuncture;
 - Needle protection against accidental puncture does not lock the needle before removing its cap;
 - Needle protection against accidental puncture after use of the needle locks the needle irreversibly.

Roman Raszewski
Director