

# Product Data Sheet

Product Name: Quicktrach I

30-04-002-1

30-04-004-1



30-04-902-1

30-04-904-1



## Order Information

REF	Description	Box
30-04-002-1	Quicktrach I Child, 2.0 mm I.D., in blister packaging, For single use, sterile	1
30-04-902-1	Quicktrach I Child, 2.0 mm I.D., in plastic tube, For single use, sterile	1
30-04-004-1	Quicktrach I Adult, 4.0 mm I.D., in blister packaging, For single use, sterile	1
30-04-904-1	Quicktrach I Adult, 4.0 mm I.D., in plastic tube, For single use, sterile	1

## Intended Use

Percutaneous cricothyrotomy set for securing the airway in case of upper airway obstruction.

## Key data

61.1.1. - 2mm

61.1.1. - 4mm

61.1.2.

REF	30-04-002-1	30-04-902-1	30-04-004-1	30-04-904-1
<b>Patient</b>	Child	Child	Adult	Adult
<b>Inner diameter needle</b>	2.0 mm	2.0 mm	4.0 mm	4.0 mm
<b>Cannula:</b> <b>Connector to patient's airway.</b> Cone in accordance to ISO 5356-1	15 mm O.D.	15 mm O.D.	15 mm O.D.	15 mm O.D.
<b>Extension Tubing:</b> <b>Connector to patient's airway.</b> Cone in accordance to ISO 5356-1	15 mm O.D.	22 mm I.D.	15 mm O.D.	22 mm I.D.
<b>Material</b> <b>Cannula</b> <b>Stopper</b> <b>Needle</b> <b>Syringe</b> <b>Necktape</b>  <b>Extension tubing</b>	ABS (acrylonitrile butadiene styrene), EVA (ethylene vinyl acetate copolymers) ASA (acrylonitrile styrene acrylate) Stainless Steel PP (polypropylene) Foam: PU (polyurethane), PA (polyamide) Hook and loop fastener: PES (polyester), PP (polyamide) Silicone                      PP                      Silicone                      PP (Polypropylene)                      (Polypropylene)			
<b>Plastic Tube dimension</b> Height Diameter	n/a	23.0 cm 4.4 cm	n/a	23.0 cm 4.4 cm
<b>Box weight</b>	140 g	140 g	145 g	145 g
<b>Box dimension</b> Length Width Height	24.5 cm 12.9 cm 3.5 cm	24.5 cm 5.7 cm 5.5 cm	24.5 cm 12.9 cm 3.5 cm	24.5 cm 5.7 cm 5.5 cm

## Product Characteristics

<b>Sterile packed:</b>	Yes <input checked="" type="checkbox"/>	No <input type="checkbox"/>	<b>RoHS compliant:</b>	Yes <input type="checkbox"/>	No <input type="checkbox"/>
Sterilization method:	EO			N/A <input checked="" type="checkbox"/>	
<hr/>					
<b>Single use product:</b>	Yes <input checked="" type="checkbox"/>	No <input type="checkbox"/>	<b>REACH SVHC included:</b>	Yes <input type="checkbox"/>	No <input checked="" type="checkbox"/>
<hr/>					
<b>Phthalates</b> requiring labelling acc. to CLP Regulation (EC) 1272/2008	Yes <input type="checkbox"/>	No <input checked="" type="checkbox"/>	<b>Use of natural rubber latex:</b>	Yes <input type="checkbox"/>	No <input checked="" type="checkbox"/>
			61.1.3. sudetyje nera latekso		
<hr/>					
<b>Shelf-life:</b>	5 years				

---

Revision: 01/12.2021

---

**VBM Medizintechnik GmbH**

Einsteinstrasse 1 | 72172 Sulz a. N. | Germany  
Tel.: +49 7454 9596-0 | Fax: +49 7454 9596-33  
e-mail: info@vbm-medical.de | www.vbm-medical.de