


# PremiCron®

## Non-absorbable sutures

### DESCRIPTION

PremiCron® is intended for use in soft tissue approximation of the wound edges to render possible an undisturbed wound healing.

PremiCron® is indicated for use in general surgery, cardiovascular and vascular indications and reconstructive and plastic surgery, when surgical practice requires the use of a non-absorbable suture material. PremiCron® is intended for adult patients regardless of gender who do not meet any contraindication.



PremiCron®

## PremiCron®

### Product overview

Polyester, braided, coated





Colour	Green or white
Origin	Synthetic
Sizes	USP 6/0 (metric 0.7) to USP 5 (metric 7)
Type of absorption	Non-absorbable
Sterilisation	Ethylene oxide / Gamma-irradiation

# PremiCron®

Needle length mm	Thread length cm colour	USP (metric)						
		0 (3.5)	1 (4)	2 (5)	3 (6)	4 (6)	5 (7)	6 (8)

1/2 circle reverse cutting needle

 	HS 37s	75 green	C0026467 13.1 p.d.	C0026468	C0026469			
	HS 40	75 green		C0026479	C0026480	B0026481		B0026482
	HS 40s	75 green			C0026390			
	HS 40ss	75 green		C0026397				
	HS 48	75 4 x 75 200 green	C0026495	C0026375 B0120070	B0026376	B0026377		B0026399 G0026203
	HS 60	75 green						B0026409

Needle length mm	Thread length cm colour	USP (metric)						
		4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)	2 (5)	3 (6)

Straight reverse cutting needle

 	GS 51	75 green		C0026505	C0026506	C0026507	C0026508	B0026509
	GS 60	75 90 green		C0026314	C0026315 C0026441	C0026266		
	GS 65	75 green			C0026287	C0026288	C0026289	







C Box of 36 units

















B Box of 24 units



















G Box of 12 units

M Box of 6 units

# PremiCron®

Needle		Thread	USP (metric)						
length mm		length cm colour	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)	2 (5)	3 (6)
1/2 circle round bodied needle									
 		75 green						C0026070	
		75 90 120 green	C0026014	C0026015  C0026002	C0026016 C0026049 C0026111	C0026017 13.2 p.d. C0026001	C0026018		
 		4 x 75 green			M0026725				
	 	4 x 75 white			M0027112				

Needle		Thread	USP (metric)							
length mm		length cm colour	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)	2 (5)	3 (6)	
1/2 circle round bodied needle										
		75	C2026814	C2026738	C2026739	C2026817				
		90 green		C2026815	C2026816					
	13.4 p.d. 13.3 p.d.									
		75			C2027839	C2027817				
		90 white			C2027816					
	 	90 Pledgets 6x3mm Firm green			C2026916					
	 	90 Pledgets 6x3mm Firm white			C2027306					
		75	C0026814	C0026738	C0026739	C0026817				
		90 green		C0026815	C0026816					
		75			C0027839	C0027817				
		90 white			C0027816					
		90 Pledgets 6x3mm Firm green		C0026915	C0026916					
		90 Pledgets 6x3mm Firm white			C0027306					
		4 x 75			M0027700	M0027009				
		8 x 75 8 x 90 green/white			M0027715 M0027713					
	 	8 x 75 8 x 90 green/white			M0027716 M0027717					
	 	4 x 75 8 x 75 8 x 90 green/white			M0027766 M0027712	M0027759				

Needle		Thread	USP (metric)							
length mm		length cm colour	6/0 (0.7)	5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)	
1/2 circle round bodied needle										
		45	C2026802	13.6 p.d. 13.5 p.d.			C2026606			
		75		C2026803	C2026804	C2026805	C2026806			
		90		C2026903	C2026904	C2026905	C2026906			
		green								
		75			C2027455		C2027457			
		white								
	 	75					C2027042			
		green								
		45						C0026606		
		75		C0026803	C0026804	C0026805	C0026806			
		90		C0026903	C0026904	C0026905	C0026906			
		green								
		75			C0027455		C0027457			
		90					C0027906			
	white									
	8 x 75						M0027775			
	4 x 90							M0027008		
green/white										
 	4 x 75						M0027764			
	8 x 75						M0027802			
green/white										
 	4 x 75						M0027088	M0027091		
	8 x 75						M0027776			
green/white										
		75					C0026257			
		green								
		8 x 75					M0027151			
		green/white								
	 	4 x 75					M0027090			
8 x 75						M0027119				
green/white										
 	8 x 75					M0027150				
	green/white									
		4 x 45				C0026740				
	take off	green								