


60.1 pd



HAEMONETICS®
THE Blood Management Company

Donor

Disposable

Therapy

REF 0980E-00

Therapeutic Plasma Exchange Closed Disposable Set, 225ml Bowl

Description

The TPE disposable set is used to collect plasma from a patient and replace it with fluid, such as Albumin, FFP, or balanced electrolyte solutions.

Replacement Solution administration can be machine controlled or optionally operator controlled through a second venous access

Technical Information

Disposable definition

Bowl

Latham centrifuge bowl, (225ml).

Needle

No needle pre-connected.

Bags

1 X Temporary substitution Bag assembly, 1000 ml

1 X Waste bag, 5 liter not pre-connected

1 X Air management bag, 600 ml

Tubing system

AC line with spike and drip chamber

Substitution line with spike and vented drip chamber.

Sterilization method

EtO


Packaging

8 units/case

Accessories

Anticoagulant solution

On the decision of the physician



Donor

Disposable

Therapy

REF 0981E-00

Therapeutic Plasma Exchange Closed Disposable Set, 125ml Bowl

Description

The TPE disposable set is used to collect plasma from a patient and replace it with fluid, such as Albumin, FFP, or balanced electrolyte solutions.

Replacement Solution administration can be machine controlled or optionally operator controlled through a second venous access

Technical Information

Disposable definition

Bowl

Latham centrifuge bowl, (125ml).

Needle

No needle pre-connected.

Bags

1 X Temporary substitution Bag assembly, 1000 ml

1 X Waste bag, 5 liter not pre-connected

1 X Air management bag, 600 ml

Tubing system

AC line with spike and drip chamber

Substitution line with spike and vented drip chamber.

Sterilization method

EtO


Packaging

8 units/case

Accessories

Anticoagulant solution

On the decision of the physician



Donor

Disposable

Therapy

REF 0971E-00

Peripheral Blood Stem Cells Closed Collection Set, 125ml Bowl

Description

The disposable set is used for the collection of peripheral blood stem cells or for therapeutic reductions of platelets or white cells.
Enhanced safety features
One-arm procedure.
Mobile (MCS+ allows bedside collection).

Technical Information

| | |
|-----------------------|-----------------------------------------------------------------------------------------|
| Disposable definition | |
| Bowl | Latham bowl (125 ml) |
| Needle | Pre-connected needle 16 G with needle guard and vacuum sampling adapter. |
| Bags | 1 X PBSC storage bag, 600 ml 1 X Plasma bag, 1000 ml (Chimney) 1 X Buffy coat bag |
| Tubing system | Anticoagulant line with spike and bacterial filter |
| Sterilization method | EtO |
| Packaging | 8 units/case |

Accessories

| | |
|------------------------|-------------------|
| Anticoagulant solution | ACD-A recommended |
|------------------------|-------------------|



60.2 pd



HAEMONETICS®
THE Blood Management Company

Donor

Disposable

Plasma

REF 0625B-00

Blow molded Centrifuge Bowl

Description

LN 0625B-00 is a centrifuge bowl used in plasmapheresis procedures.


Technical Information

| | |
|-----------------------|---------------|
| Disposable definition | |
| Bowl Capacity | 275 ml |
| Sterilization method | Radiation |
| Packaging | 25 units/case |

Accessories

| | |
|------------------------------|------------------------------------------------------------------------|
| Plasma collection disposable | 00690-00, 00692-00, 00699-00, SC690-00, SC692-00, PB1000SPV, PB1000LLV |
| Donor harness | 00620-00, 0620E-00 |
| Related disposables | 00698-00, 03018-01 |
| Anticoagulant | Sodium citrate 4%, (ratio 1:16) |

60.2 pd



HAEMONETICS®
THE Blood Management Company

Donor

Disposable

Plasma

REF 00620-00

Plasma Collection Donor Harness

Description

LN 00620-00 is a plasma collection donor harness


Technical Information

| | |
|-----------------------|------------------------------------|
| Disposable definition | |
| Harness length | 183 cm (72 inch) |
| Bowl connector | For connection to centrifuge bowl. |
| Needle connection | Luer connector |
| Tubing system | AC line with spike |
| Sterilization method | EtO |
| Packaging | 100 units/case |

Accessories

| | |
|-----------------|--------------------|
| Centrifuge bowl | 0625B-00, 625HS-00 |
|-----------------|--------------------|

66.2 pd



HAEMONETICS®
THE Blood Management Company

Donor

Disposable

Plasma

REF 00690-00

Plasma Collection Bag

Description

LN 00690-00 is a plasma collection disposable for the collection of plasma for fractionation.

Technical Information

| | |
|-----------------------|------------------------------------|
| Disposable definition | |
| Capacity | 1000 ml |
| Bowl connector | For connection to centrifuge bowl. |
| Sterilization method | EtO |
| Packaging | 102 units/case |

Accessories

| | |
|-----------------|---------------------|
| Centrifuge bowl | 00625-00B, 625HS-00 |
| Donor harness | 00620-00, 0620E-00 |