

# Ordering Codes

## Small Diameter Reload with Tri-Staple™ 2.0 Technology

The chiped 8 mm Signia™ Small Diameter Reload (with 6.3 mm Tip) can be used with any of our manual and powered handles. Pair it with our Signia™ Stapler for improved access and smart features like:

- Visual and audible feedback before and during firing.
- The ability to adjust firing speeds to optimize staple formation in variable tissue conditions.

Order-Code	Color	Description	Box Qty
<b>SIGSDS30CTV</b>		<b>Small Diameter Curved Tip Reload, 30 mm Vascular Tissue</b> Shaft Length 15 cm, Staple Heights 2.0 mm, 2.0 mm	<b>6</b>
<b>SIGSDS30CTVT</b>	<input type="checkbox"/>	<b>Small Diameter Curved Tip Reload, 30 mm Vascular / Thin Tissue</b> Shaft Length 15 cm, Staple Heights 2.5 mm, 2.5 mm	<b>6</b>
<b>SIGSDL45CTVT</b>	<input type="checkbox"/>	<b>Small Diameter Curved Tip Reload, 45 mm Vascular / Thin Tissue</b> Shaft Length 24 cm, Staple Heights 2.5 mm, 2.5 mm	<b>6</b>

3.6, 3.7 Kasetės darbinės dalies ilgis 15 cm. Atidarytos kabutės aukštis 2 mm