

Invitrogen™

## Barrier (Filter) Tips, 20 µL size

Catalog number: AM12645

Related applications: [RNA Extraction](#) [PCR](#) [Real Time PCR \(qPCR\)](#)



**Catalog Number**

**AM12645**

**Unit Size**

10 racks

**Price (USD)**

[Contact Us](#) ▶



Save to list

[Shop now](#)

Find the Essentials lab products you need

Ambion® Barrier tips prevent cross-contamination of samples by means of an aerosol filter. They are RNase- and DNase-free and are ideal for handling RNA, PCR applications, and for working with radioisotopes or toxic chemicals. Ten racks containing 8 x 12 tips each are provided. They are compatible with Eppendorf®, Finnpiquette, Oxford Benchmate®, Pipetman™, Proline®, Rainin, Socorex, and other brands of pipettors. Each lot of tips is subjected to rigorous testing for RNase contamination and is certified to be nuclease-free.

For Research Use Only. Not for use in diagnostic procedures.

## Specifications

<b>Volume (Metric)</b>	20 µL
<b>Filtered</b>	Filtered
<b>Sterility</b>	Sterile
<b>Autoclavable</b>	Autoclavable
<b>Format</b>	Rack
<b>Tip Style</b>	Universal Fit
<b>Packaging</b>	10 x 96 tips
<b>Shipping Condition</b>	Room Temperature
<b>For Use With (Equipment)</b>	Universal Fit
<b>Product Line</b>	Invitrogen

## Contents & Storage

Store at room temperature.

[Shop now](#)

Find the Essentials lab products you need

## Certificates

Search by lot number or partial lot number

Lot #	Certificate Type	Date	Catalog Number(s)
 <a href="#">30121295</a>	Certificates of Analysis	Feb 22, 2023	AM12645
 <a href="#">19321402</a>	Certificates of Analysis	Feb 02, 2023	AM12645
 <a href="#">10621402</a>	Certificates of Analysis	Aug 03, 2021	AM12645
 <a href="#">10221402</a>	Certificates of Analysis	Jul 21, 2021	AM12645
 <a href="#">10421402</a>	Certificates of Analysis	Jul 10, 2021	AM12645

5 results displayed, search above for a specific certificate

[Request a Certificate](#)

## Product Information

[Product Information](#)

[Shop now](#)

Find the Essentials lab products you need

---

## Manuals



[RNase-free Tips](#)