

Pump Protection Filter /Heidelberger Extension Line 200 cm



Sets offer a simple, cost-effective patient connection method, specifically designed for use with the latest digital Oocyte Aspiration Pumps.

HYDROPHOBIC BARRIER provides pump protection by preventing accidental aspiration of fluid into the device.

REMOVES CONTAMINATION RISK from irregular or uncontrolled decontamination of multiple use consumable filters.

LIGHTWEIGHT PATIENT TUBE SET: 2.5m lightweight patient tube set with integral luer lock connector to ensure secure, reliable connection to all standard luer-lock fitting aspiration needle sets.

DESIGNED SPECIFICALLY to match the performance characteristics of current digital aspiration pumps.

Filter Set for Oocyte Aspiration Pump 0.2µm filter set, 2.5m patient connection tube with male luer connector. Supplied sterile, for single use in cartons of 10 units