

RESTORATIVE ACCESSORIES
FINISHING & POLISHING
FINISHING & POLISHING

ceram.x® gloss

finishing & polishing system

Product Description

ceram.x gloss finishing & polishing system is a single use product, to finish and polish the surface of composite and compomer restorations in a two-step technique.

A new generation of polishing instruments for a perfect finalization of your ceram.x universal restoration.

- Exceptionally intuitive: one shape for all tooth structures
- Excellent gloss: especially harmonized system with diamond abrasives for outstanding polishing results
- Extremely fast: high gloss within just a few minutes

ceram.x gloss Finishing and Polishing

Refill

30 x Finishing
30 x High-Gloss Polishing

60701612

ceram.x gloss & ceram.x universal

Combi Pack

52 x 0.25 g ceram.x universal Compules Tips
30 x ceram.x gloss Finishing
30 x ceram.x gloss High-Gloss Polishing

Shade A2 (=Vita® shades A2, B2, D2)
Shade A3 (=Vita® shades A3, C2, D3, D4)

60701613
60701614

Combi Pack

4 x ceram.x universal Syringes (2 x A2, 2 x A3)
30 x ceram.x gloss Finishing
30 x ceram.x gloss High-Gloss Polishing

60701615

