




## INJECTION NEEDLE



Primary packaging: Pouch  
Secondary packaging: cardboard box

### Description:

vienkartinio naudojimo (pažym ta simboliu

- Triple bevel needle for an easy penetration,
- Ergonomic handle to control needle movements,
- Transparent catheter for liquid flow visualization.

tinka endoskopo kanalui: Ø2 mm

| References       | Description | Catheter      |             | Needle        |             | Handle color | Minimum delivery quantity |
|------------------|-------------|---------------|-------------|---------------|-------------|--------------|---------------------------|
|                  |             | Diameter (mm) | Length (cm) | Diameter (mm) | Length (mm) |              |                           |
| CJ-A3-23-230-21* | Novaject    | 2.3           | 230         | 21 / 0.8      | 3           | ●            | 10                        |
| CJ-A3-23-230-23* |             | 2.3           | 230         | 23 / 0.6      | 3           | ●            |                           |
| CJ-A3-23-230-25* |             | 2.3           | 230         | 25 / 0.5      | 3           | ●            |                           |
| CJ-A-18-165-25   |             | 1.8           | 165         | 25 / 0.5      | 4           | ●            |                           |
| CJ-A-23-160-19   |             | 2.3           | 160         | 19 / 1.1      | 4           | ●            |                           |
| CJ-A-23-230-19   |             | 2.3           | 230         | 19 / 1.1      | 4           | ●            |                           |
| CJ-A-23-230-21   |             | 2.3           | 230         | 21 / 0.8      | 4           | ●            |                           |
| CJ-A-23-230-23   |             | 2.3           | 230         | 23 / 0.6      | 4           | ●            |                           |
| CJ-A-23-230-25   |             | 2.3           | 230         | 25 / 0.5      | 4           | ●            |                           |
| CJ-A5-23-230-21* |             | 2.3           | 230         | 21 / 0.8      | 5           | ●            |                           |
| CJ-A5-23-230-23  |             | 2.3           | 230         | 23 / 0.6      | 5           | ●            |                           |
| CJ-A5-23-230-25  |             | 2.3           | 230         | 25 / 0.5      | 5           | ●            |                           |

\* Special request, delivery time on demand

CE Mark  
0197




23G 5mm

### NovaJect

REF CJ-A5-23-230-23

≥2.8mm

Model: VDK-IN-23-230-2305-A




#### Disposable Injection Needle

FR Aiguille d'injection / EN Disposable Injection Needle / DE Einmal-Injektionsnadel / ES Aguja de inyección desechable / NL Injectienaald voor eenmalig gebruik / PT Agulha de Injecção Descartável / IT Ago iniezione monouso / DK Engangs injektionsnål / GR Βελόνα εγχυστικής χρήσης / SE Engangs Injektionsnål


- ant pakuot s pažym tas produkto galiojimo laikas m nesiais;

LOT 2022092007VIN




REF CJ-A5-23-230-23  
LOT 2022092007VIN

2025-08



REF CJ-A5-23-230-23  
LOT 2022092007VIN

2022-09



REF CJ-A5-23-230-23  
LOT 2022092007VIN

Jiangsu Vitekang Medical Science & Technology Co., Ltd.  
No. 52 Guosheng Road  
Wujin Economic Development Zone Changzhou  
213149 Jiangsu, China

Prolinx GmbH  
Brehmstr. 56, 40239 Duesseldorf,  
Germany

REF CJ-A5-23-230-23  
LOT 2022092007VIN

VYTIL  
161 avenue Gallieni  
93170 Bagnolet  
info@vytil.com



Version 3 - Date 2021-03-11

(01)06970057576039(17)250831(10)2022092007VIN