

Optilene®

Non-absorbable sutures


10 p.d.

DESCRIPTION

Optilene® is a sterile non-absorbable monofilament surgical suture made of polypropylene and polyethylene.

Available in blue for a better visibility.

Optilene® is offered with RacePack® and LongPack® packaging presentations to reduce the memory effect. Other special presentations, including PremiPatch® PTFE pledgets, are also available.



Optilene®

Product overview

Polypropylene and polyethylene, monofilament, uncoated

| | |
|----------------------|---|
| Colour | Blue |
| Origin | Synthetic |
| Sizes | USP 10/0 (0.2 metric) to USP 2 (5 metric) |
| Types of degradation | Non-absorbable |
| Sterilisation | Ethylene oxide |

| Needle length mm | Thread length cm colour | USP (metric) | | | | | | |
|--|---|--------------|--------------|------------|------------|------------|----------|----------|
| | | 5/0 (1) | 4/0 (1.5) | 3/0 (2) | 2/0 (3) | 0 (3.5) | 1 (4) | 2 (5) |
| 1/2 circle round bodied needle | | | | | | | | |
|  HR 30 EasySlide |  | 45 | | | | | | |
| | | 75 | | | C3090046 | C3090069 | | |
| | | 90 | | | | C3090047 | C3090048 | C3090049 |
| | | 120 | | | | C3095927 | C3090477 | |
| | | blue | | | C3090144 | | | |
|  2xHR 30 EasySlide |  | 90 | | | | | | |
| | | 120 | | | C3090001 | C3090002 | | |
| | | blue | | | C3090462 | | | |
|  HR 30s EasySlide |  | 75 | | | | | | |
| | | 100 | | | | C3090092 | C3090093 | |
| | | blue | | | | C3090096 | | |
|  HR 30ss EasySlide |  | 75 | | | | | | |
| | | blue | | | | C3090094 | | |
|  HR 37 EasySlide |  | 75 | | | | | | |
| | | 90 | | | C3090050 | C3090051 | C3090052 | C3090053 |
| | | blue | | | | C3095789 | | |
|  2xHR 37 EasySlide |  | 90 | | | | | | |
| | | 120 | | | C3090942 | C3090943 | C3095054 | |
| | blue | | | C3090981 | C3090973 | 10.1 p.d. | | |
| | | | | | | 10.2 p.d. | | |
|  HR 37s EasySlide |  | 2 x 90 | | | | | | |
| | | blue | | | G3095727 | | | |
|  HR 37s EasySlide |  | 75 | | | | | | |
| | | 120 | | | C3090150 | C3090151 | C3090152 | C3090153 |
| | | blue | | | | C3090154 | | C3095340 |
|  2xHR 37s EasySlide |  | 90 | | | | | | |
| | | 120 | | | C3090982 | | | |
| | | blue | | | C3090971 | C3090972 | | |

| Needle length mm | Thread length cm colour | USP (metric) | | | | | | | |
|--|---|--|------------------|--------------|------------|--------------|------------|------------|----------|
| | | 8/0 (0.4) | 7/0 (0.5) | 6/0 (0.7) | 5/0 (1) | 4/0 (1.5) | 3/0 (2) | 2/0 (3) | |
| 1/2 circle round bodied needle | | | | | | | | | |
|  2xHR 22 | | 60 | | | | C3090084 | | | |
| | | 75 | | | | | C3090863 | C3090086 | |
| | | 90 | | | | C3090975 | C3090976 | C3090860 | C3090978 |
| |  | 120 | | | | | C3090977 | C3090978 | |
| | | blue | | | | | C3090876 | | |
|  2xHR 22b |  |  Pledgets 6x3mm Firm | 90 | | | | C3090979 | | |
| |  |  Pledgets 6x3mm Soft | 90 | | | | C3095336 | | |
| |  | Long pack | 2 x 75 2 x 90 | | | | G3095620 | G3095616 | |
|  2xHR 22s |  | | 75 | | | | | C3095508 | |
| |  | | 90 | | | | C3095512 | C3095513 | |
|  2xHR 22fb |  | | 90 | | | C3096035 | | | |
| | | | blue | | | | | | |
|  2xHR 26ff |  | | 90 | | | | C3095541 | | |
| | | | blue | | | | | | |
|  2xHR 26f |  | | 75 | | | | C3096004 | C3096003 | C3096021 |
| |  | | 90 | | | | C3095611 | C3096017 | C3096020 |
| |  | | 120 | | | | C3096030 | C3096028 | C3096026 |
|  2xHR 26fb |  |  Pledgets 6x3mm Firm | 90 | | | | C3096910 | | |
| | | | blue | | | | | | |
|  HR 26 |  | | 90 | | | | C3095606 | C3095596 | |
| | | | blue | | | | | | |
|  HR 26 |  | | 75 | | | C3090039 | C3090040 | C3090042 | |
| | | | blue | | | | | C3090041 | |

10.3 p.d.

Optilene®

| Needle length mm | Thread length cm colour | USP (metric) | | | | | | | |
|--|---|--|---|---------------------|------------|------------|----------|----------|--|
| | | 5/0 (1) | 4/0 (1.5) | 3/0 (2) | 2/0 (3) | 0 (3.5) | 1 (4) | 2 (5) | |
| 1/2 circle round bodied needle | | | | | | | | | |
|   |  | 75 | | | | | | C3090043 | |
| | | 90 | | C3090140 | | | | | |
| | | 120 blue | | | | | C3090444 | C3095712 | |
|   |  | 75 | C3090906 | C3090707 | C3090708 | C3090087 | | | |
| | | 90 | C3090905 | C3090907 | C3090908 | C3090909 | | | |
| | | 120 blue | | C3095709 | C3090709 | C3095710 | | | |
| | | | | 10.5 p.d. 10.4 p.d. | | | | | |
| | | |  | 90 | | C3090910 | C3090838 | | |
|   |  | 90 | | C3095347 | | | | | |
| | | |  | 90 | | C3095347 | | | |
| | | |  | 2 x 90 blue | | G3090633 | G3090631 | | |
|   |  | 75 | | C3095721 | | C3095723 | | | |
| | | 90 blue | | C3095557 | C3095558 | | | | |
|   |  | 2 x 90 blue | | G3095993 | | | | | |
| | |   |  | 90 | | C3095647 | | | |
| 120 | | | | C3095647 | C3095635 | | | | |
| blue | | | | | | | | | |
|   |  | 75 | | | C3095722 | | | | |
| | | 90 blue | | C3095688 | | | | | |
|   |  | 75 | | | C3090015 | C3095442 | C3095443 | C3095452 | |
| | | 90 blue | | | | | | C3095458 | |
|   |  | 120 blue | | | C3090959 | C3090956 | | | |
|   |  | 75 90 blue | | | | C3090032 | C3095444 | C3090449 | |

| | | | |
|-------------------|-------------------|-------------------|------------------|
| C Box of 36 units | B Box of 24 units | G Box of 12 units | M Box of 6 units |
|-------------------|-------------------|-------------------|------------------|

Optilene®

| Needle length mm | Thread length cm colour | USP (metric) | | | | | | | |
|---|---|--------------|--------------|--------------|------------|--------------|------------|------------|----------|
| | | 8/0 (0.4) | 7/0 (0.5) | 6/0 (0.7) | 5/0 (1) | 4/0 (1.5) | 3/0 (2) | 2/0 (3) | |
| 1/2 circle round bodied needle | | | | | | | | | |
|  |  | 75 | | | | C3090501 | C3090503 | C3090502 | |
| | | 90 | | | C3090898 | C3090901 | C3090902 | C3090903 | C3090900 |
| | 120 | | | | | C3090904 | | | |
| | blue | | | | | 10.8 p.d. | 10.6 p.d. | | |
| |  | 90 | | | | | C3090837 | | |
|  | 90 | | | | | C3095327 | | | |
|  |  | 2 x 75 | | | | | | G3095183 | |
| | | 2 x 90 | | | | G3090642 | G3090194 | G3090644 | |
| | | 4 x 90 | | | | G3095868 | | | |
| | | blue | | | | | | | |
|  |  | 75 | | | | C3095460 | C3095461 | | |
| | | 90 | | | | C3095492 | C3095493 | C3095494 | |
| | | 120 | | | | | | C3095994 | |
| | | blue | | | | | | | |
|  |  | 90 | | | | C3095660 | C3095643 | | |
| | | blue | | | | | | | |
|  |  | 90 | | | C3096157 | C3096156 | C3090729 | | |
| | | 120 | | | | | | | C3096155 |
| | | blue | | | | | | | |
|  |  | 90 | | | | C3090223 | C3090338 | C3090426 | |
| | | blue | | | | | | | |
|  |  | 75 | | | | | C3095503 | | |
| | | 90 | | | | | C3095526 | | |
| | | blue | | | | | | | |
|  |  | 75 | | | | | C3096009 | C3096011 | C3096019 |
| | | 90 | | | | | C3096010 | C3096008 | C3096018 |
| | | 120 | | | | | C3096029 | C3096027 | C3096025 |
| | blue | | | | | | | | |
|  |  | 90 | | | | | C3096114 | | |
| | | blue | | | | | | | |
|  |  | 60 | | | | | | C3095720 | |
| | | 75 | | | | | C3090029 | C3090030 | |
| | | blue | | | | | | | |

| | | | |
|-------------------|-------------------|-------------------|------------------|
| C Box of 36 units | B Box of 24 units | G Box of 12 units | M Box of 6 units |
|-------------------|-------------------|-------------------|------------------|

Optilene®

| Needle length mm | Thread length cm colour | USP (metric) | | | | | | | |
|--|--|--------------------------|--------------|--------------|----------------------|--|--|----------------------------------|----------------------------------|
| | | 8/0 (0.4) | 7/0 (0.5) | 6/0 (0.7) | 5/0 (1) | 4/0 (1.5) | 3/0 (2) | 2/0 (3) | |
| 3/8 circle round bodied needle | | | | | | | | | |
|   |  | 45 | | | | | | | |
| | | 75 blue | | | |  10.7 p.d. | C3096031 C3090120 | C3095057 | |
|  |  | 60 | | | | | | | |
| | | 75 90 120 blue | | | | C3095918 C3090719 C3090919 | C3090720 C3090920 | C3095920 C3095921 | C3095092 |
|   |   | 4 x 75 2 x 90 blue | | | | | | | |
| | | | | | G3090394 G3095659 | G3090659 | | | |
|   |  | 90 120 blue | | | | | | | |
| | | | | | C3095242 | C3095656 | | | |
|   |  | 90 blue | | | | | | | |
| | | | | | | C3095683 | | | |
|   |  | 60 blue | | | | | | | |
| | | | | | | C3096154 | | | |
|   |  | 75 | | | | | | | |
| | | 90 | | | | | | | |
| | | 120 | | | | C3090821 | C3090722 C3090822 C3095112 C3096130 | C3090723 C3090823 C3095420 | C3090721 C3090824 C3095114 |
| | | 75 blue | | | | | | | |
|   |  | 90 blue | | | | | | | |
| | | | | | | | C3096153 | | |
|   |  | 75 blue | | | | | | C3090132 | |
| | | | | | | | | | |
|   |  | 75 blue | | | | | | C3090724 C3090725 | |
| | | | | | | | | | |
| 3/8 circle round bodied needle | | | | | | | | | |
|   |  | 75 blue | | | | | | C3090082 C3090083 | |
| | | | | | | | | | |

| | | | |
|-------------------|-------------------|-------------------|------------------|
| C Box of 36 units | B Box of 24 units | G Box of 12 units | M Box of 6 units |
|-------------------|-------------------|-------------------|------------------|

| Needle length mm | Thread length cm colour | USP (metric) | | | | | | |
|---|---|--------------|--------------|--------------|------------|--------------|------------|------------|
| | | 8/0 (0.4) | 7/0 (0.5) | 6/0 (0.7) | 5/0 (1) | 4/0 (1.5) | 3/0 (2) | 2/0 (3) |
| 1/2 circle round bodied needle | | | | | | | | |
|  2xHR 13f |  | 45 | | | C3097120 | | | |
| | | 75 | | | C3097123 | | | |
| | | 75 | | | | | C3095185 | |
| | | 90 | | | | C3095792 | | |
| | | blue | | | | | | |
|  | 2 x 75 | | | G3095127 | G3095128 | | | |
| | 4 x 75 | | | | G3095126 | | | |
|  | 75 | | | C3095666 | C3095653 | | | |
| | blue | | | | | | | |
| 1/2 circle round bodied needle | | | | | | | | |
|  new |  | 75 | | | C3097060 | | | |
| | | blue | | | | | | |
| 1/2 circle round bodied needle | | | | | | | | |
|  HR 13 |  | 45 | | | | C3096033 | | |
| | | 75 | | | C3090499 | C3090500 | | |
|  2xHR 13 |  | 45 | | | C3090853 | C3090839 | | |
| | | 60 | | | C3095935 | C3090934 | C3090933 | |
| | | 75 | | C3090060 | C3090953 | C3090954 | | C3090045 |
| | | 90 | | | | 10.9 p.d. | C3090750 | |
| | | blue | | | | | | |
|  | 2 x 60 | | | | G3090146 | | | |
| | 2 x 75 | | | G3090625 | | G3095010 | | |
|  2xHR 13b |  | 75 | | | C3095448 | C3095449 | | |
| | | blue | | | | | | |
|  2xHR 17f |  | 75 | | | | C3095501 | | |
| | | 90 | | | | C3095992 | C3096001 | |
| | | blue | | | | | | |
|  | 2 x 90 | | | | | G3095521 | | |
| | blue | | | | | | | |
|  HR 17 |  | 45 | | | | | C3096032 | C3090038 |
| | | 75 | | | | C3090012 | C3090013 | C3090014 |
| | | blue | | | | | | |
|  HR 17 |  loop | 90 | | | | | C3090413 | |
| | | blue | | | | | | |

Optilene®

| Needle length mm | Thread length cm colour | USP (metric) | | | | | | |
|---|---|--------------|----------------------|--------------|---|--------------|------------|------------|
| | | 8/0 (0.4) | 7/0 (0.5) | 6/0 (0.7) | 5/0 (1) | 4/0 (1.5) | 3/0 (2) | 2/0 (3) |
| 1/2 circle round bodied needle | | | | | | | | |
|  new HR 10  |  75 blue | | C3097401 | C3097402 | | | | |
|  new 2xHR 10  |  60 75 blue | | C3097488 C3097867 | C3097868 | | | | |
| |  2 x 75 blue | | | G3097878 | | | | |
|  new 2xHR 10b  |  75 blue | | | C3097007 | | | | |
| 1/2 circle round bodied needle | | | | | | | | |
|  HR 10  |  75 blue | | C3090401 | C3090402 | | | | |
|  2xHR 10  |  60 75 blue | | C3090488 C3090867 | C3090868 | C3090866 C3090869 10.10 p.d. | | | |
| |   Pledgets 3x3mm Firm | 60 blue | | | | C3095950 | | |
| |  2 x 60 2 x 75 blue | | | | G3095878 | G3090770 | | |
|  2xHR 10b  |  75 blue | | | C3096007 | | | | |
| 1/2 circle round bodied needle | | | | | | | | |
|  new 2xHR 13ff  |  75 blue | | | C3097341 | C3097346 | | | |
| |   Pledgets 3x3mm Soft | 75 blue | | | | C3097006 | | |
| |  75 blue | | | C3095341 | C3095346 | | | |
| |   Pledgets 3x3mm Soft | 75 blue | | | | C3096006 | | |

| | | | |
|-------------------|-------------------|-------------------|------------------|
| C Box of 36 units | B Box of 24 units | G Box of 12 units | M Box of 6 units |
|-------------------|-------------------|-------------------|------------------|

Optilene®

| Needle length mm | Thread length cm colour | USP (metric) | | | | | | |
|---|---|--------------|--------------|--------------|------------|--------------|------------|------------|
| | | 8/0 (0.4) | 7/0 (0.5) | 6/0 (0.7) | 5/0 (1) | 4/0 (1.5) | 3/0 (2) | 2/0 (3) |
| 3/8 circle round bodied needle | | | | | | | | |
|  2xDR 10b  |  | 60 | | C3096131 | | | | |
| | | blue | | | | | | |
| 3/8 circle round bodied needle | | | | | | | | |
|  DR 13  |  | 75 | | | C3098930 | C3097931 | | |
| | | blue | | | | | | |
|  2xDR 13  |  | 45 | | C3098885 | C3098886 | C3098887 | | |
| | | 60 | | C3097894 | C3097895 | C3097896 | | |
| | | 75 | | C3097905 | C3097906 | C3097907 | | |
| |  | 90 | | | C3097938 | C3097939 | | |
| | | blue | | | | | | |
| | | 75 | | | | | C3098173 | |
|  2xDR 13b  |  | 2 x 45 | | | G3097957 | | | |
| | | 2 x 60 | | | G3097970 | | | |
| | | 4 x 60 | | | G3097571 | | | |
| | | 75 | | | G3097472 | | | |
| | | 2 x 75 | | | G3097979 | G3098980 | | |
| | | 4 x 75 | | | G3097582 | | | |
|  DR 13  |  | 60 | | | C3097926 | C3097609 | | |
| | | 75 | | | C3097191 | C3097547 | | |
| | | 90 | | | C3097839 | C3097840 | | |
| | | blue | | | | | | |
| 3/8 circle round bodied needle | | | | | | | | |
|  DR 13  |  | 75 | | | C3095930 | C3095931 | C3095932 | |
| | | blue | | | | | | |

| | | | |
|-------------------|-------------------|-------------------|------------------|
| C Box of 36 units | B Box of 24 units | G Box of 12 units | M Box of 6 units |
|-------------------|-------------------|-------------------|------------------|

| Needle length mm | Thread length cm colour | USP (metric) | | | | | | |
|---|---|--------------|--------------|-----------------------|------------|--------------|------------|------------|
| | | 8/0 (0.4) | 7/0 (0.5) | 6/0 (0.7) | 5/0 (1) | 4/0 (1.5) | 3/0 (2) | 2/0 (3) |
| 3/8 circle round bodied needle | | | | | | | | |
|  |  | 45 | | C3097614 | C3097615 | | | |
| | | 60 | | C3097714 | C3097715 | C3097522 | | |
| | | 75 | | C3097914 | C3097915 | C3097913 | | |
| | | 90 | | C3097799 | C3097815 | | | |
| | | blue | | 10.13 p.d. 10.12 p.d. | | | | |
| |  | 75 | | | C3097306 | | | |
| | | 2 x 45 | | | G3097636 | | | |
| | | 4 x 45 | | G3097369 | | | | |
| | | 2 x 60 | | G3097611 | G3097613 | | | |
| | | 4 x 60 | | G3097371 | | | | |
|  | 2 x 75 | | G3097810 | G3097811 | | | | |
| | 4 x 75 | | G3097373 | | | | | |
| | blue | | | | | | | |
| | | | | | | | | |
| 3/8 circle round bodied needle | | | | | | | | |
|  |  | 45 | | C3097164 | | | | |
| | | 60 | | C3097062 | C3097816 | | | |
| | | 75 | | C3097812 | C3097843 | | | |
| | | 90 | | | | | | |
| | | blue | | 10.14 p.d. | | | | |
| |  | 75 | | | C3095306 | | | |
| | | 2 x 45 | | | G3090636 | | | |
| | | 4 x 45 | | G3090369 | | | | |
| | | 2 x 60 | | G3090611 | G3090613 | | | |
| | | 4 x 60 | | G3090371 | | | | |
|  | 2 x 75 | | G3095810 | G3095811 | | | | |
| | 4 x 75 | | G3090373 | | | | | |
| | blue | | | | | | | |
| | | | | | | | | |
|  |  | 45 | | C3090164 | | | | |
| | | 60 | | C3090062 | C3095816 | | | |
| | | 75 | | C3095812 | C3095843 | | | |
| | | blue | | | | | | |
|  |  | 60 | | C3096262 | | | | |
| | | 75 | | | C3095679 | | | |
| | | blue | | | | | | |

| Needle length mm | Thread length cm colour | USP (metric) | | | | | | | |
|--|---|---|--------------|--------------|------------|--------------|------------|------------|--|
| | | 8/0 (0.4) | 7/0 (0.5) | 6/0 (0.7) | 5/0 (1) | 4/0 (1.5) | 3/0 (2) | 2/0 (3) | |
| 3/8 circle round bodied needle | | | | | | | | | |
| 2xDR 8fb  new  |  | 45 | C3097824 | C3097826 | | | | | |
| | | 60 | C3097829 | C3097832 | | | | | |
| 3/8 circle round bodied needle | | | | | | | | | |
| 2xDR 8f  |  | 45 | C3090107 | C3090079 | | | | | |
| | | 60 | C3090105 | C3090088 | | | | | |
| | |  | 2 x 45 | G3090118 | G3090124 | | | | |
| | | | 2 x 60 | G3090121 | G3090127 | | | | |
| | | | 2 x 75 | | G3090130 | | | | |
| 3/8 circle round bodied needle | | | | | | | | | |
| 2xDR 8fb  |  | 45 | C3095824 | C3095826 | | | | | |
| | | 60 | C3095829 | C3095832 | | | | | |
| 3/8 circle round bodied needle | | | | | | | | | |
| 2xDR 8  |  | 12 | | C3097870 | | | | | |
| | | 45 | C3097703 | C3097685 | C3097686 | C3097684 | | | |
| | | 60 | C3097887 | C3097888 | C3097889 | | | | |
| | | 75 | C3097004 | C3097885 | C3097886 | | | | |
| | |  | 2 x 45 | G3097361 | G3097362 | G3097399 | | | |
| | | | 2 x 60 | G3097363 | G3097364 | | | | |
| | | | 75 | | G3097089 | | | | |
| | | | 2 x 75 | G3097892 | G3097080 | | | | |
| | | | blue | | | | | | |
| | | 3/8 circle round bodied needle | | | | | | | |
| 2xDR 8b  new |  | 45 | C3097805 | | | | | | |
| | | 60 | C3097817 | C3097819 | C3097821 | | | | |
| 3/8 circle round bodied needle | | | | | | | | | |
| 2xDR 8  |  | 12 | | C3095870 | | | | | |
| | | 45 | C3090703 | C3090685 | C3090686 | C3090684 | | | |
| | | 60 | C3090887 | C3090888 | C3090889 | | | | |
| | | 75 | C3090004 | C3090885 | C3090886 | | | | |
| | |  | 2 x 45 | G3090361 | G3090362 | G3090399 | | | |
| | | | 2 x 60 | G3090363 | G3090364 | | | | |
| | | | 75 | | G3090089 | | | | |
| | | | 2 x 75 | G3090892 | G3095080 | | | | |
| | | | blue | | | | | | |

10.15 p.d.