


The disposable cable is available to order in multiples of 1 each. The Model 4051L is the alligator style cable that is compatible with all leads, including the IS4/DF4 with the testing sleeve.

The testing sleeve is available in Optisure leads packaging and is planned to be including in other IS4 leads when globally approved. The IS4/DF4 testing sleeve (Model EX3151) is available through customer service as an accessory in multiples of 1 each.



The 4161 is another option of disposable cables but has a hook instead of an alligator connection.



**Model 4161**

Please contact Technical Service Support for clarification on the PSA cable options that are available

- **Phone:** +46 8 474 4147
- **Email:** [technical.support@sjm.com](mailto:technical.support@sjm.com)

# SJM™ Peel-Away Introducer

14 cm Sheath

Introducer Kit

5 F – 16 F

## Product Highlights

- Proprietary materials improve insertion characteristics and reduce vessel trauma
- Close tolerance extrusion and proprietary tipping process facilitate tracking on a guidewire
- Reliable peeling characteristics
- Sheath can be totally occluded without kinking to prevent air aspiration
- Di-Lock™ feature secures dilator in sheath during insertion

## Ordering Information

SJM™ Peel-Away Sheath, Di-Lock™ Dilator, 12 cc syringe, 18 ga. XTW Needle and 50 cm Guidewire with 3 mm “J” (5 units per box).

REORDER NUMBER	FRENCH SIZE	MAXIMUM GUIDEWIRE DIAMETER (IN)	SHEATH USABLE LENGTH (CM)
405104	6	.038	14
405108	7	.038	14
405112	8	.038	14
405116	9	.038	14
405118	9.5	.038	14
405120	10	.038	14
405122	10.5	.038	14
405124	11	.038	14
405128	12	.038	14
405136	14	.038	14

