



### IOL Delivery Systems


The MONARCH® IOL Delivery System combines a reusable titanium handpiece and a sterile single use cartridge for enhanced implantation of the AcrySof® IOL. Advanced design enables the surgeon to view and verify lens orientation throughout the implantation process. Simplified loading, controlled consistent delivery, and ease of implantation are the benefits that MONARCH® IOL Delivery System can provide.



**MONARCH® II IOL Injector** 8065-977771



**MONARCH® III IOL Injector** 8065-977773



**INTREPID® AUTOSERT® Handpiece** 8065-751755

### MONARCH® Cartridges

10/box



50	<b>A Cartridge</b>	8065-977757
49	<b>B Cartridge</b>	8065-977758
47	<b>C Cartridge</b>	8065-977762
48	<b>D Cartridge</b>	8065-977763

### MONARCH® II Loading Forceps 560.01

The ALCON®/GRIESHABER® MONARCH® II Loading Forceps is for fully controlled handling of the AcrySof® Single-Piece and AcrySof® Natural IOLs from its packaging into the MONARCH® Delivery System.

[ 3 ]