

ROEKO Aspirator tips

The plastic aspirator tips for single use with an opening of 4.8 mm are ideal for moisture evacuation during all standard procedures where sterility is not a requirement. Supplied with an autoclavable Double Adapter to fit on all equipment.



Aspirator tip without green tip	
Tip-opening \varnothing 6.5 mm / End-opening \varnothing 4.8 mm	
100 pcs*	461 048

*includes 1 Double Adapter

ROEKO Concorde II

The smooth contours of the Concorde II aspirator tips mean they are comfortable for the patient and easy to clean. Concorde II is made from high quality, durable plastic.

- low noise suction
- flange retracts tongue and cheek gently
- matt surface prevents light reflection
- can be sterilized in autoclave or chemiclave up to 134 °C
- the tips are available with \varnothing 11 mm and \varnothing 16 mm for all standard equipment. The 11 mm tips are supplied with an autoclavable Adapter.



Size		
\varnothing 11 mm*	10 pcs	450 211
\varnothing 16 mm	10 pcs	450 216

*includes 1 Adapter

Adapter and Double Adapter

To ensure that ROEKO aspirator tips can be used on all standard aspirating equipment, Adapters are available. All Adapters are autoclavable and chemiclavable up to 134 °C.



Double Adapter	463 000
The Double Adapter can be used to fit Surgitip, Surgitip-micro and the non sterile aspirator tips onto \varnothing 11 mm and \varnothing 16 mm equipment.	

Adapter		
For connecting Concorde II and Surgitip canulas to 11 and 16 mm aspirating equipment		
Adapter	\varnothing 6.5 > \varnothing 11 mm	463 025
Adapter	\varnothing 11 > \varnothing 16 mm	463 011