

## FLOW COMPOSITE

**Master-Dent® Flow Composite** is a light cured flowable hybrid resin restorative recommended for Class III, Class IV and Class V restorations. Combines the Bis-GMA chemistry with a non-slumping flowability for easier placement and contouring. Flows into the preparation versus traditional packing and manipulation. Dispensed through a syringe with an application tip for increased accuracy and fingertip control.



| REF              |   |             |                                    |
|------------------|---|-------------|------------------------------------|
| <b>Intro Kit</b> |   |             |                                    |
| 19-414-AS        | 4 x 2Gm syringes (1 ea. A1, A2, A3, C3) w/ tips |             |                                    |
| REF              |   | REF         |                                    |
| Refill Kits      | 4 x 2Gm syringes w/ Tips                        | Refill Kits | 4 x 2Gm syringes w/ Tips           |
| 19-414-A1        | Shade A1 w/ tips                                | 19-414-C2   | Shade C2 w/ tips                   |
| 19-414-A2        | Shade A2 w/ tips                                | 19-414-C3   | Shade C3 w/ tips                   |
| 19-414-A3        | Shade A3 w/ tips                                | 19-414-C4   | Shade C4 w/ tips                   |
| 19-414-A3.5      | Shade A3.5 w/ tips                              | 19-414-D1   | Shade D1 w/ tips                   |
| 19-414-A4        | Shade A4 w/ tips                                | 19-414-D2   | Shade D2 w/ tips                   |
| 19-414-B1        | Shade B1 w/ tips                                | 19-414-D3   | Shade D3 w/ tips                   |
| 19-414-B2        | Shade B2 w/ tips                                | 19-414-G    | Shade Gingival w/ tips             |
| 19-414-B3        | Shade B3 w/ tips                                | 19-414-O    | Shade Opaque w/ tips               |
| 19-414-C1        | Shade C1 w/ tips                                | 19-411      | Flow Composite refill tips 100/bag |

## SPECTRON™ TRA MICROHYBRID LIGHT CURE UNIVERSAL COMPOSITE

Spectron-TRA™ is a light-cured, resin-based radiopaque all purpose restorative that is a highly polishable, 76% filled micro-hybrid composite with an average filler particle size of less than 7 microns. The material applies easily and is ideally suited for anterior and posterior applications.



| REF   |   |                 |          |
|---|---|-----------------|----------|
| <b>Standard Kit</b>   |   |                 |          |
| 19-350  | 4 x 4.5 Gm Syringes (1 ea. A1, A2, A3, B2),<br>7 Gm Ultimate Bond, 3 mL Etch Syringe, Accessories |                 |          |
| <b>Trial Size Kit</b>   |   |                 |          |
| 19-370  | 2 x 4.5 Gm Syringes (1 ea. A2, A3),<br>7 Gm Ultimate Bond, 3 mL Etch Syringe, Accessories         |                 |          |
| REF   |   | REF             |          |
| 4.5 Gm Syringes   |   | 4.5 Gm Syringes |          |
| 19-351  | Shade A1  | 19-361          | Shade C1 |
| 19-352  | Shade A2  | 19-362          | Shade C2 |
| 19-360  | Shade A3  | 19-363          | Shade C3 |
| 19-353  | Shade A3.5  | 19-364          | Shade C4 |
| 19-355  | Shade A4  | 19-369          | Shade D1 |
| 19-354  | Shade B1  | 19-371          | Shade D2 |
| 19-357  | Shade B2  | 19-367          | Shade D3 |
| 19-368  | Shade B3  | 19-372          | Shade D4 |
| 19-373  | Shade B4  |                 |          |
| 19-215  | Ultimate Bond 7 Gm  |                 |          |
| * <b>COMPULE REFILLS:</b> 20 Compules @0.25 Gm ea. Same REF numbers as above (syringe refills) but add a C for compules such as #19-351-C |   |                 |          |