





















# Aesculap® YASARGIL Aneurysm Clip System

## Phynox Standard Clips

PHYNOX STANDARD

PERMANENT						PERMANENT					
Art. No. Art. Nr.		Blade length Maullänge mm	Maximal opening maximale Öffnung mm	Closing force Schließkraft N g		Art. No. Art. Nr.		Blade length Maullänge mm	Maximal opening maximale Öffnung mm	Closing force Schließkraft N g	
FE740K		7.0	6.2	1.47	150	FE762K		10.2	7.5	1.77	180
FE750K		9.0	7.0	1.77	180	FE782K		13.7	8.7	1.96	200
FE760K		11.0	7.8	1.77	180	FE786K		15.3	9.8	1.96	200
FE780K		15.0	9.2	1.96	200	FE744K		5.4	5.5	1.47	150
FE792K		17.5	10.6	1.96	200	FE754K		6.4	6.1	1.77	180
FE790K		20.0	11.4	1.96	200	FE764K		8.0	6.5	1.77	180
4.3.3. FE840K		25.0	13.3	1.77	180	FE784K		10.3	7.4	1.96	200
FE742K		6.5	6.0	1.47	150	FE823K		7.4	5.4	1.77	180
FE752K		8.3	6.8	1.77	180	FE824K		8.4	5.8	1.77	180