



### 36p.d Savaime išsipleiantis tulžies latak stentas

Sign in [Forgot Details?](#)

Email Address

.....




Remember Me

[Register FREE](#)

PPE | [Interventional Radiology](#) | [Imaging](#) | [Cardio-thoracic](#) | [Orthopaedics](#) | [Anaesthesia](#)

[Browse By Manufacturer](#) > [Cook Medical](#) > [Zilver Stent](#)



## Zilver Stent

Added Sep 6, 2007

Editor Rating

Manufactured by [Cook Medical](#)

User Rating

[Add your Rating](#)

Search by [Advanced search](#)

Enter your search term

[Suggest a Device](#)

[Add to My Medical Devices](#)

[My Medical Devices](#)

Categorised under:  
[Interventional Radiology](#) > [Vascular](#) > [Stents \(non covered\)](#)

[Description](#) | [Reviews \(1\)](#) | [Comments \(0\)](#) | [Videos \(0\)](#) | [Evidence \(0\)](#) | [News/Editorial \(1\)](#) | [Tips & Tricks \(0\)](#) | [Alerts \(0\)](#)

### Description

[From Which Medical Device](#)

Used as an adjunct to percutaneous transluminal angioplasty (PTA) in the treatment of symptomatic vascular disease of the iliac arteries.

- Self-expandable stent made of nitinol.
- Flexible, slotted tube which provides strength and flexibility upon deployment.
- Supplied with multiple gold markers on each end for precise placement in the lumen.
- Pre-loaded in a 5.0 Fr OD Flexor® introducer - designed to provide maximum flexibility without kinking or compressing.
- Recommended 5.0 Fr sheath/7.0 Fr guiding catheter size.
- Accepts .018 inch diameter wire guide.

Supplied sterile in peel-open packages. Intended for one-time use.

### Have you used this device?

Your opinion matters to others - [rate this device](#) or [add a comment](#)

[Similar Devices](#) 1-5 of 25

Sort By: [Newest/Recently Added](#) ▾