

LigaV® Ligating Clips

Product Description Card

Ref. no.: 0301-01S 0301-01M 0301-01ML 0301-01L 0301-01ML04 0301-01S20 0301-01M20


Grena Ltd, 1000 Great West Road, Brentford,
Middlesex, TW8 9 HH, United Kingdom

Contact information:
Phone: +44 131 208 207 6
Fax: +44 131 777 80 77


 0197

Product name:	LigaV® Ligating Clips	Product generic group	Ligating Clips
Class	IIb – rule 8	GMDN Code	35649

Intended use: Marking and/or ligating of any tissue structures or vessels during an operation for haemostasis or marking purposes. The product can be used only by the qualified medical personnel.

Product description: LigaV® ligating clips are implantable medical devices made of medical grade titanium, which ensure proper level of biocompatibility. The LigaV® clip is formed into a V-shape from a rectangular wire with the appropriate serrations on the inner side of the clip that enhance its grip on the tissue and enable its proper nourishment. The outer serration helps to grip the clip in the applicator's jaws and gives better stability of the clip.

Product technical drawing



Product dimensions*					
REF:	Clips per cartridge	Clip shape in the cartridge	Opened clip length (outer) after loading into applier	Opened clip width (internal) after loading into applier	Closed clip length
0301-01S	6	V shaped	3,20 mm	2,10 mm	3,70 mm
0301-01M	6	V shaped	4,70 mm	2,90 mm	5,40 mm
0301-01ML	6	V shaped	8,00 mm	5,20 mm	9,00 mm
0301-01L	6	V shaped	11,90 mm	5,40 mm	12,30 mm

* the presented dimensions are average

Material specification:	
Clip	Titanium
Cartridge for clips S, M	Acrylonitrile butadiene styrene
Cartridge for clips ML, L	Nylon
Special materials content:	
DEHP	NO
Latex	NO

Certificates:	
EC certificate MDD 93/42/EEC Annex II	
EN ISO 13485:2012/AC 2012	
ISO 13485:2003 under CMDCAS	
EN ISO 9001:2008	
Notified body	TÜV Rheinland

Additional information:	
Recommended tissue thickness	
S	0,30-1,50 mm
M	1,00-2,50 mm
ML	2,50-4,00 mm
L	3,50-7,50 mm

Packing Specification*	
Primary packaging:	
Packing type	Blister
Packing material	Tyvek/PET
Quantity	1 cartridge per blister
Secondary packaging:	
1. Packing material	Carton box
Size S, M (6 pcs.)	50 blisters per box
Size ML, L (6 pcs.)	20 blisters per box
Size ML (4 pcs.)	20 blisters per box
2. Packing material	Cardboard box
Size S, M (6 pcs.)	600 blisters per box
Size ML, L (6 pcs.)	240 blisters per box
Size ML (4 pcs.)	240 blisters per box

* not applicable for 01S20 and 01M20

Sterilization:	
Sterile?	YES
Sterilized by	EO
Re-sterilization	NO – Disposable use only
Shelf life	5 years

Symbol	Description	Symbol	Description	Symbol	Description	Symbol	Description
	Caution, consult accompanying documents		Keep dry		Consult instructions for use		Manufacturer
	Do not reuse		Do not use if package is damaged		Do not re-sterilize		Sterilized by EtO