

6 BPJ - Embedding

1| Made of polyoxymethylene, a high quality technical plastic, these embedding cassettes are totally resistant to solvents and decalcifying agents used in histological techniques


2| Their special design allows a maximum flow through and proper drainage of all solutions used in the process of dehydration and paraffin embedding of the tissue samples.

3| Lids are supplied attached to the cassette and can be opened and closed as often as necessary; they always relock securely without danger of specimen loss.

4| Each cassette presents 3 large labelling surfaces for proper sample identification.

Atsparios ksileno poveikiui

6 Universal cassettes with lid




1| Frontal labelling area 45° inclined

Code	Colour	Pore diameter	Pack
BPJ001	White	1x5 mm	250
BPJ002*	Blue	1x5 mm	250
BPJ003*	Pink	1x5 mm	250
BPJ004*	Yellow	1x5 mm	250
BPJ005*	Green	1x5 mm	250

*Codes availables until end of stock

6 Biopsy cassettes with lid



1| Frontal labelling area 45° inclined

Code	Colour	Pore diameter	Pack
BPJ006	White	1x1 mm	250
BPJ007*	Blue	1x1 mm	250
BPJ008*	Pink	1x1 mm	250
BPJ009*	Yellow	1x1 mm	250
BPJ010*	Green	1x1 mm	250

*Codes availables until end of stock

6 CCassettes without lid



Audinių paruošimas automatinuose audinių paruošimo procesoriuose ir audinių parafininio bloko paruošimas

Pritaikytos žymėti paprastu pieštuku

1| Frontal labelling area 45° inclined

2| The base allow an optimal adherence to the paraffin block.

3| Used preferably for making paraffin blocks, they fit securely to any microtome chuck adapter.

Code	Colour	Pore diameter	Pack
BPJ011	White	2 mm	500
BPJ012*	Blue	2 mm	500
BPJ013*	Pink	2 mm	500
BPJ014*	Yellow	2 mm	500
BPJ015*	Green	2 mm	500

*Codes availables until end of stock

Tinkamos parafininių pjūvių ruošimui pusiau automatiniais mikrotomais



6 BPJ - Embedding

6 Stainless lid for cassetes

- 1 | Stainless steel
- 2 | For embedding cassetes without lid.

Code	Description
BPJ039	Stainless Steel lid for cassettes

