

Erika Ciegieniė | Sentios

From: Maite de Andrés <export4@auxilab.es>
Sent: Friday, April 21, 2023 3:46 PM
To: Erika Ciegieniė | Sentios
Subject: Embedding cassettes BPJ011

Flag Status: Flagged

Hello Erika,
We inform the ref. BPJ011 matches with these specifications.

Best regards,

Maite de Andrés (Mrs)
Export Sales Assistant

Tel: + (34)948310513 - EXT. 227

export4@auxilab.es

www.auxilab.es



De: Erika Ciegieniė | Sentios <erika.ciegieniė@sentios.lt>

Enviado el: jueves, 20 de abril de 2023 10:19

Para: Maite de Andrés <export4@auxilab.es>

Asunto: Embedding cassettes BPJ011

Hello,

Can you confirm that Embedding cassettes (without lid, B/500, white), Ref. BPJ011 comply all the following requirements:

- Field of application: Histology.

- Purpose: Tissue preparation in automatic tissue preparation processors and tissue paraffin block preparation.
- Suitable for preparation of paraffin sections with semi-automatic microtomes.
- Resistant to xylene.
- Suitable for marking with a simple pencil.

Thank you.

Best regards,
Erika Čiegienė
Manager

UAB „Sentios“

Address: A. Juozapaviciaus pr. 104b, Kaunas, LT-45108, Lithuania

Email: erika.ciegiene@sentios.lt

Phone: +370 670 11421

www.sentios.lt

sentios