

When your job is making the right decisions

Product Range





multiFiltrate acute dialysis therapy system

The range of therapies offered by the multiFiltrate acute dialysis machine are shown below

- CVWH** Continuous veno-venous haemofiltration
- CVVHD** Continuous veno-venous haemodialysis
- CVVHDF** Continuous veno-venous haemodiafiltration
- HV-CVWH** High volume continuous veno-venous haemofiltration
- SCUF** Slow continuous ultrafiltration
- MPS** Membrane plasma separation
- HP** Haemoperfusion

multiFiltrate with full set of lines



Heater bags allow quick installation and deinstallation

System shown fully set up

multiFiltrate cassette

Dialysate line

Substitution line

Heparin syringe

The setup of the tubing systems is clearly shown and explained in steps on the screen of the multiFiltrate



Waste bags

Type	Description	Connection	Units per box	Art. No.
Waste bag 2000 mL	Waste bag with female Luer-Lock connector		100	501 509 1
Filtrate bag 10 L	Filtrate collection bag with draining cock, male Luer-Lock connector		40	502 901 1
Filtrate bag 10 L Single use	Filtrate collection bag with Luer-Lock male, single use		40	502 903 1
PF waste bag 5000 mL	Waste bag with male Luer-Lock connector for single use		50	501 510 1

Accessories

Type	Description	Connection	Units per box	Art. No.
Pressure monitor line	Pressure monitor line with filter, connector male Luer-Lock, length of 30 cm	 	100	501 463 1
Pressure measurement line 60 cm	Pressure monitor line with filter, connector male Luer-Lock, length of 60 cm, for arterial pressure measurement	 	100	501 915 1
Check valve	Check valve Luer-Lock male/ Luer-Lock female	 	10	850 080 1
Scissors clamp	For clamping of the lines		5	284 024 1
Heparin syringe 50 mL	Fresenius Injectomat syringe, Luer-Lock, latex-free		50	900 071 1
Heparin syringe 50 mL	Fresenius P-syringe, latex-free, Luer-Lock		50	900 076 1
Dialyser exchange set	Tubing system for exchange of dialyser during treatment, short-circuit AV-set, 2 lines for priming of dialyser		75	501 800 1
Freka-Flex Infusion set	Infusion set with roller clamp and drip chamber	 	150	288 901 1

Saline solutions

Type	Presentation	Units per box	Units per pallet	Art. No.
Isotonic saline solutions (1 litre contains sodium chloride 9.0 g, Na ⁺ 154 mmol/L, Cl ⁻ 154 mmol/L)	30 × 250 mL		1,440 Units	131 067 1
	15 × 500 mL		720 Units	131 865 1
	10 × 1000 mL		480 Units	131 866 1

If products are categorised as pharmaceuticals, availability depending on pharmaceutical product registration.

Fresenius Medical Care Deutschland GmbH · 61346 Bad Homburg

**Vilnius University Hospital
Santaros Clinics**

**Fresenius Medical Care
Deutschland GmbH**
Product Unit MFTpro & Genius
Critical Care & Ventures
Care Enablement
Else-Kröner-Straße 1
61352 Bad Homburg
Germany
www.freseniusmedicalcare.com

August 07th, 2023

**Subject: Tender no. 680686 Products for hemofiltration procedures
and other disposable medical aids (6136)**

To whom it may concern

We, Fresenius Medical Care Deutschland GmbH, located in Else-Kroener-Strasse 1, 61352 Bad Homburg, Germany, the official sales and distribution entity of Fresenius Medical Care AG & Co. KGaA, located in Else-Kroener-Strasse 1, 61352 Bad Homburg, Germany, herewith confirm that our product **multiEffluent Bag 10 L (art. no. 5029011)** is a **single use medical device** intended to be used at the acute dialysis machines multiFiltratePRO and multiFiltrate.

Best regards



I.V. Dirk Etdorf
Director



I.V. Vera Hering-Dorn
Senior Product Manager

KOPIJA TIKRA 