

Surgical Loop

For organ retraction and more

Description

Surgical Loop are retraction tapes made of radiopaque silicone, cotton or polyester. The surgical tapes are used for the intraoperative isolation, marking, mechanical support and looping of organs, blood vessels, tendons and nerves, or for tying off the umbilical cord of the new born infant.

Indications

Surgical Loop is used to:

- Retract nerves, tendons, arteries, veins, ureters, vessels and other structures to access the operation site
- Differentiate and/or identify organs or areas while operating
- Clamp

Advantages

Security

- Strong and reliable
- Single-use
- Atraumatic
- Silicone radiopaque Surgical Loop



Convenience

- Extremely smooth surface
- Pliable
- No tissue adherence
- Easy, convenient and safe packaging presentation
- Distinct colours to ease differentiation
- Various product dimensions to suit all purposes

Surgical Loop

Product overview

Structure	Tube, braided tape
Colour	White, blue, green, yellow, red, violet
Chemical composition	Silicone, cotton, polyester
Coating	Uncoated
Origin	Synthetic, natural
Sizes	1.2 - 2.5 mm (silicone) 4 - 8 mm (cotton) 3.2 - 4 mm (polyester)
Type of absorption	Non-absorbable
Sterilisation	Ethylene oxide

Surgical Loop

For organ retraction and more

Radiopaque silicone tube

Diameter	Length	Units / box	Colour / code			
						
1.5 mm	2 x 45 cm	24	B1095048	B1095021	B1095030	B1095013
1.5 mm	75 cm	24	B1095242	B1095226	B1095234	B1095218
2.5 mm	2 x 45 cm	24	B1095145	B1095129	B1095137	B1095110
2.5 mm	75 cm	24	B1095544	B1095528	B1095536	B1095510
1.2 mm	2 x 45 cm	24	B1095552	The thread is attached to a DRN 25 needle		

Cotton braided tape **2 p.d.**

Width	Length	Units / box	Colour / code			
						
4 mm	2 x 38 cm	24	B1094920			
4 mm	75 cm	24	B1095641	B1095633	B1095625	B1095650
8 mm	75 cm	24	B1095676		B1095668	
4 mm	75 cm	24	B1095617	The tape is attached to a GS 65 needle		
				GS 65 ▼ 		

Polyester braided tape

Width	Length	Units / box	Colour / Code
			
3.2 mm	2 x 75 cm	24	B1094904
4 mm	75 cm	24	B1094912

C Box of 36 units	B Box of 24 units	G Box of 12 units	M Box of 6 units
-------------------	-------------------	-------------------	------------------

Surgical Loop **Skirtas organų atitraukimui**

Surgical Loop tai juostelė pagaminta iš silikono, medvilnės ir poliesterio. Chirurginė juostelė skirta intraoperacinei izoliacijai, markiravimui, mechaniniam palaikymui ir kilpai ant organų, kraujagyslių, sausgyslių ir nervų, ir naujagimių bambagyslei perrišti.

Chirurginės atitraukimo juostelės pagamintos iš silikono, medvilnės ir poliesterio

- Atskirti organų vietą operacijos metu ir atitraukti ar prilaikyti nervams, sausgyslėms, arterijoms, venoms, ureteriui ar kraujagyslėms ir suteikti lengva priėjimą operacijos metu
- Užspaudimui
- Užrišti naujagimių bambagyslę.

Surgical Loop garantuoja ir numato:

- Lankstumą ir atsparumą
- Lygų paviršių be audinių sukibimo
- Patrauklus spalvų pristatymas lengvam atskyrimui
- Skirtingas produktų dydis pasirinktinai nuo rinkinio
- Silikoninė medžiaga garantuoja paciento saugumą.

Sudėtis – silikonas, medvilnė, poliesteris

Struktūra – vamzdelio ar pinto tipo

Dengimas – nedengtas

Spalva – raudota, mėlyna, geltona, balta, žalia ir violetinė

Storis – Vamzdelių – 1,2mm, 1,5mm ir 2,5mm

Juostelių – 3,2mm, 4mm ir 8 mm

Rezorbicija – nesirezorbuoja

Indikacijos – organų atitraukimui

Kontraindikacijos – nežinomos

Sterilizacija – etileno oksidas