


Easy to use  
Locks in large quantities of fluid without dripping

# DryMax Floor mat 2.4



DryMax Floor mat 2.4 is a handy sized superabsorbent floor mat for busy operating rooms. It can be placed under the operating table, in front of or beside the feet of the surgical team.

DryMax Floor mat 2.4 absorbs and retains a large amount of liquid: up to 3.5 litres of saline solution or 9 litres of tap water. During absorption, the mat swells and its contents turn to gel.


The mat is dual-sided and can be placed on floors that are already wet. If required, additional mats can be added to provide extra absorption.

Proper fluid management promotes Health and Safety as it reduces the risk of slips and trips and provides a cleaner and safer working environment.

As DryMax Floor mat 2.4 locks in the contaminated fluids without dripping, it helps reduce the risk of spreading contaminants. In this way, the mat contributes to hygiene control and staff safety.

Cleaning times and costs are also significantly reduced as fluids are absorbed and retained by the mat.

Item	REF	Pcs/ p.d.1	Product size	Absorbency		Packaging
				Tap water	NaCl 0,9%	
DryMax Floor mat 2.4	F20241	100	72 cm x 37 cm	9 l	3,5 l	25 pcs in a plastic bag x 4 = 100 pcs/box
DryMax Floor mat 2.4 bulk	F20241-B	100	72 cm x 37 cm	9 l	3,5 l	100 pcs



**DISPOSAL**  
Treat as clinical waste, combustion forms carbon dioxide and water



Single-use

Absorbs water-based liquids such as saline, urine, amniotic fluids and blood. In the case of liquids containing bodily fluids and proteins, the absorption capacity will decrease compared to the saline data provided.

# OPERATION DryMax<sup>®</sup> Floor mat 2.4



- Operating Rooms
- Examination Rooms
- Gynaecology
- Urology
- Emergency rooms
- Ambulances
- Body bags
- Staff Safety
- Hygiene Control
- Faster turnover time
- Easy to use

**ABSORBEST**

Manufactured by Absorbent AB Sweden  
www.absorbent.se | info@absorbent.se