



## Multifunctional chair JORDAN T3e (1.)



### Technical parameters:

1. For examination of patients, if necessary, the chair is transformed into a lying position **(2.)**
2. It consists of 3 sections: **(3.)**
  1. Head - back, buttocks, calves. **(3.1.)**
  2. The frame of the chair is made of steel, painted with polyester powder paint. **(3.2.)**
3. Head - back section inclination angle from 0° to 80° **(4.)**
4. Angle of inclination of the calf part from 0° to - 70° **(5.)**
5. Seat section inclination angle from 0° to - 15° **(6.)**
6. Chair length (horizontal position) 1900 mm **(7.)**
7. Seat width (excluding armrests) 700 mm **(8.)**
8. Maximum safe permissible load 200 kg **(9.)**

9. Chair control: **(10.)**
1. Leg part. With the help of a gas spring **(10.1.)**
  2. The back part. With the help of an electric motor **(10.2.)**
  3. Sitting part. With the help of an electric motor **(10.3.)**
  4. Trendelenburg position. **(10.4.)**
10. Adjustable armrests (angle rotation, height change) **(11.)**
11. Chair chassis - with 4 wheels. wheel diameter 75 mm, 2 wheels with brakes. **(12.)**



RESI Třeboň spol. s r.o. , Novohradská 1153, 379 01 Třeboň, Czech Republic  
Tel.: +420 384 724 292, [www.resi.cz](http://www.resi.cz)

