

# Polypectomy



**US** endoscopy

*listening...and delivering solutions\**

# Polypectomy

US Endoscopy offers an extensive line of polypectomy devices. From a purpose-built cold snare to a monofilament EMR device, we offer a product solution to accommodate a variety of needs.

## **Exacto® cold snare**

Studies have shown that using the purpose-built Exacto® cold snare results in a higher complete resection rate compared to traditional braided wire snares<sup>1</sup>.

*wire diameter is 33% thinner than traditional braided wires allowing for full circumferential resection of polyps*

*shield shape and mini snare size offers more control for placement and resection of small polyps*

*reduces polyp "fly away" from the resection site making it easier to collect specimen for pathology<sup>2</sup>*



Placement of the Exacto® snare around polyp



Cold snaring with the Exacto® snare



Clean cut at resection site with Exacto® snare

<sup>1</sup> A. Horiuchi, K. Hosoi, M. Kajiyama, N. Tanaka, K. Sano, D. Graham. "Prospective, randomized comparison of 2 methods of cold snare polypectomy for small colorectal polyps." GIE Journal Online, 2015.  
<sup>2</sup> S. Din, A. Bal, S. Riley, P. Kitsanta, S. Johal. "Cold snare polypectomy: Does snare type influence outcomes?" Digestive



## Lariat® snare

The multi-functional design of the Lariat® snare helps address a variety of resection needs including piecemeal resection and multiple polyps.

*three distinct sizes and shapes eliminate the need to exchange snares*

*snare loop design maintains its shape and integrity after multiple resections*



## iSnare® system

A unique, multi-functional system ideal for endoscopic mucosal resection (EMR) and the removal of polyps.

*inject or snare at any time throughout the procedure without device exchanges while maintaining scope position*

*facilitates immediate injections in response to bleb absorption/dissipation or bleeding*

*available with oval or hexagonal snare (snares also available stand-alone in 15 and 25mm)*



# **NEW!** **Histolock™ resection device**

The Histolock™ resection device features a stiff monofilament wire to help recruit and resect flat tissue during complex polypectomy procedures.

*twisted wire that lays flat on the mucosa during resection*

*mini snare size that offers easy maneuverability for controlled placement and precise resection*

**Histolock**  
*resection device*

## **NEW!** **Histolock™ resection kits**

The Histolock™ resection device is also available in two convenient EMR kit options which include the needed polypectomy and tissue retrieval devices for challenging cases...



5-piece kit

- *Histolock™ resection device*
- *Carr-Locke injection needle*
- *Roth Net® retriever - polyp\**
- *eTrap® polyp trap*
- *Mio® medical device organizer*

*\*(included in BX00711873 kit only)*

# Polyp Retrieval

## ***Roth Net® retrievers***

From capturing large polyps to securing multiple polyp fragments, the Roth Net® retriever is a device you can count on for consistent performance. Our all-encompassing line of retrieval nets are designed to...

*quickly and easily retrieve multiple, large or piecemeal polyps*

*preserve specimens for histology*

*allow for continued examination with clear visibility*



Roth Net® retrieval of large polyp



Roth Net® retrieval of multiple polyps/polyp fragments



EMR specimen removal with Roth Net® retriever



## eTrap<sup>®</sup> polyp trap

Allows for quick and accurate retrieval of multiple specimens, while safeguarding clinicians from unnecessary exposure to biomaterials.



*collect multiple polyps with uninterrupted suction*

*simple transport and preservation of specimens*

*magnifying window allows easy visualization of the collected specimen*

*includes 2 trays*




## Reveal<sup>™</sup> distal attachment cap

The Reveal<sup>™</sup> distal attachment cap is designed to facilitate the endoscope's view during polypectomy procedures.

*helps avoid 'red out' by maintaining a constant distance between the endoscope and mucosa*

*can be used with other polypectomy devices, including snares and needles*



side hole allows for constant drainage of fluid



# product information

product number	description	needle gauge/ projection	sheath diameter (mm)	length (cm)	approximate snare size (mm)	units/ box
<b>polypectomy snares</b>						
BX00711018	ShortThrow™ snare - mini	-	2.4	230	12.5 x 30	10
BX00711019	ShortThrow™ snare - standard oval	-	2.4	230	25 x 54	10
BX00711085	iSnare® system - oval	23 gauge x 5mm	3.0	230	25 x 40	5
BX00711086	iSnare® system - oval	25 gauge x 5mm	3.0	230	25 x 40	5
BX00711088	iSnare® system - hexagonal	25 gauge x 5mm	3.0	230	25 x 40	5
BX00711089	iSnare® system - hexagonal	23 gauge x 5mm	3.0	230	25 x 40	5
BX00711101	ShortThrow™ snare - standard oval	-	2.4	350	25 x 54	5
BX00711113	Traxtion® snare	-	2.4	230	25 x 40	10
BX00711115	Exacto® cold snare	-	2.4	230	9.0	10
BX00711117	Histolock™ resection device	-	2.4	230	14 x 27	10
BX00711871*	Histolock™ resection device	-	2.4	230	14 x 27	3
BX00711872	Histolock™ resection kit-4 piece	-	-	-	-	3
BX00711873	Histolock™ resection kit-5 piece	-	-	-	-	3
BX00711119	Lariat® snare	-	2.4	230	oval: 30x50 hex: 10x28 diamond: 6x15	10
BX00711120	Rotator™ snare - standard	-	2.4	230	25 x 54	10
BX00711121	Rotator™ snare - mini	-	2.4	230	12.5 x 30	10
BX00711770	Reveal™ distal attachment cap (sterile)	8.9-9.5mm endoscope tip OD	11.35mm cap OD	4mm	-	10
BX00711771	Reveal™ distal attachment cap (sterile)	9.5-9.9mm endoscope tip OD	11.8mm cap OD	4mm	-	10
BX00711772	Reveal™ distal attachment cap (sterile)	9.8-10.5mm endoscope tip OD	12.4mm cap OD	4mm	-	10
BX00711773	Reveal™ distal attachment cap (sterile)	10.9-11.8mm endoscope tip OD	13.4mm cap OD	4mm	-	10
BX00711774	Reveal™ distal attachment cap (sterile)	12.6-13.2mm endoscope tip OD	15.0mm cap OD	4mm	-	10
BX00711775	Reveal™ distal attachment cap (sterile)	13.4-13.9mm endoscope tip OD	15.7mm cap OD	4mm	-	10
<b>polypectomy retrieval</b>						
BX00711050	Roth Net® retriever	-	2.5	230	30 x 60	5
BX00711052	Roth Net® retriever	-	2.5	230	30 x 60	10
BX00711099	eTrap® polyp trap	-	-	-	-	25
BX00711150	Roth Net® retriever - polyp	-	2.5	230	30 x 60	5
BX00711151	Roth Net® retriever - polyp	-	2.5	230	30 x 60	10
BX00711158	Roth Net® retriever - mini polyp	-	2.8	230	25 x 45	10
BX00711162*	Roth Net® retriever	-	2.5	230	30 x 60	10
BX00711165*	Roth Net® retriever - maxi	-	3.0	230	40 x 80	5
BX00711168*	Roth Net® retriever - mini	-	1.8	160	20 x 45	5
BX00715050/ BX00715060*	Roth Net® Platinum® retriever - universal	-	2.5	230	40 x 55	5
BX00715150/ BX00715160*	Roth Net® Platinum® retriever - polyp	-	3.0	230	40 x 55	5

\*sterile/hot currently for sale in the U.S.

**NEW!**  
**NEW!**  
**NEW!**  
**NEW!**

**NEW!**



US Endoscopy is a subsidiary of STERIS Corporation.  
© 2015 US Endoscopy or its affiliates. All rights reserved.  
Listening...and delivering solutions®, the US Endoscopy listening...and delivering solutions design, and all marks denoted with ® or ™  
are registered with the U.S. Patent and Trademark Office, or trademarks owned by US Endoscopy.

5976 Heisley Road  
Mentor, OH 44060  
phone (1) 440 / 639-4494  
fax (1) 440 / 639-4495  
customer service (1) 800 / 769-8226

[www.usendoscopy.com](http://www.usendoscopy.com)

761130 Rev. A