



[Back to Catalog](#)

View My List 0

LinkSeq™ Wipe Test Kit (compatible with Applied BioSystems™ Standard Block or Bio-Rad CFX96 Touch™)

Catalog #
9850R-SR

Product Type
Ancillary Reagents

Add to List



Details

This product is for research use only

Regulatory: RUO: For Research Use Only. Not for use in diagnostic procedures.

Technique: Real-time PCR-SSP

Use: Used to detect genomic and amplicon contamination. Compatible with Applied Biosystems™, Standard Block or Bio-Rad CFX96 Touch™.

Storage Temperature: -35° to -15° C

Package Includes: 12 Strips, Buffer, Polymerase, Enzyme, Control DNA

Product Support Files

Search

Lot/Version
Type
English
Clear All

| Title | Type | Lot/Version | Publish Date | Download |
|---|-------------------------|-------------|--------------|----------|
| LinkSeq Wipe Test Kit, Lot K4162, Certificate of Analysis | Certificate of Analysis | K4162 | 09/17/2021 | |
| LinkSeq 96 Fast Well UDG, Amp and Melt Template for use with QuantStudio 6 Software Version 1.3 | Software Analysis Files | | 11/12/2021 | |
| LinkSeq Wipe Test Kit, Lot K4257, Certificate of Analysis | Certificate of Analysis | K4257 | 07/13/2022 | |
| LinkSeq Wipe Test Kit, Lot K4372, Certificate of Analysis | Certificate of Analysis | K4372 | 04/11/2023 | |
| LinkSeq Wipe Test Kit, 2 Samples/Strip, K3999, Certificate of Analysis, RUO | Certificate of Analysis | K3999 | 12/28/2020 | |
| LinkSeq 96 Fast Well Melt Only Template for use with QuantStudio 6 Software Version 1.3 | Software Analysis Files | | 02/11/2021 | |
| LinkSeq 96 Fast Well Amp and Melt Template for use with QuantStudio 6 Software Version 1.3 | Software Analysis Files | | 02/11/2021 | |

Kiekvienas rinkinys supakuotas atskirai