



# DATA SHEET

REF: TBB600 – TRANSFER BLOOD BAG



## Product Description

Transfer Blood Bag



## Tubing Characteristics

Inner Diameter: 0.30 cm  
Outer Diameter: 0.45 cm  
Thickness: 0.15 cm

## Sterilization

Ethylene Oxide

## Packaging

1. Individual Package: PET/PP compounded vacuum bag

- |                                    |                             |
|------------------------------------|-----------------------------|
| Qty per vacuum bag:                | 1 unit                      |
| 2. Shipping Package:               | Carton                      |
| Qty per carton:                    | 100 units (100 vacuum bags) |
| 3. Expiry Date:                    | 3 years                     |
| 4. Recommended storage conditions: |                             |

Preserve in carton box protected from light and store in a cool place, store between 0°C to 30°C.

**Note:** Statements contained in this document are for reference only. For contractual information please refer to products labels.

## GENERAL MANUFACTURING PROCEDURE IN COMPLIANCE WITH:

ISO 13485, EC Directive 93/42/EEC Annex II Article 3

Symbols:

	Red Blood cells container		Plasma container
	Buffy Coat Container		Sampling site
	Whole blood container		Non-pyrogenic fluid path
	Platelets container		Film plasticizer DEHP
	Caution		Store at 0°C -30°C
	Consult instructions for use		Expiry date
	Single use only		Sterilized by Ethylene Oxide
	Do not use if package is damaged		Lot number
	Do not re-sterilize		Catalogue Number
	Manufacturer		Manufacturing date
	CE mark and Identity No. of Notifying body. The product complies with the requirements of the Council Medical Products 93/42/EEC.		