



EDAN INSTRUMENTS INC.	
No.15 Jinhui Rd., Jinsha Community, Kengzi Subdistrict, Pingshan District, Shenzhen, 518122 P.R.China	
FAX :0086-755-26898330	Email: info@edan.com.cn www.edan.com
EDAN\QR\IM25\A	
CX12 Patient Monitor Rev 3.3.3 Year:2025	
CX12 Patient Monitor	Standard Configuration including



Main Unit
 13.3 inch, Rotary Knob
 3/5-lead ECG, RESP, NIBP, Edan SpO2, 2-TEMP, PR
 Arrhythmia, OxyCRG Display, Large Font Display, Trend Display
 240-Hour Trend Data, 48 hours Full Disclosure Waveform
Optional Configurations:
 Full-touch screen, 6/12-lead ECG, 2-IBP, CO2 (EDAN G2 Sidestream), C.O., Thermal Recorder, Built-in Wi-Fi, HDMI & Multi-function interface, Bi-directional Communications with Central Monitoring System

PART NO.	DESCRIPTION	Quantity	
83.60.263157	Standards		
Standard accessories:			Kiekis
01.57.471330	NIBP cuff, EDAN, Adult, Reusable, 27-35cm (has to work with EDAN NIBP tube) NIBP manžetė, EDAN, suaugusiems, vidutinio dydžio, daugkartinio naudojimo, 27-35 cm (turi veikti su EDAN NIBP prijungimo žarnele)	1 vnt	12.3.
01.57.471329	NIBP Cuff, Small Adult, 20.5cm-28cm, reusable, with EDAN logo (has to work with EDAN NIBP tube) NIBP manžetė, EDAN, suaugusiems, maža, daugkartinio naudojimo, 20.5cm -28cm (turi veikti su EDAN NIBP prijungimo žarnele)	1 vnt.	12.3.
02.57.225029	SpO2 Finger Sensor, Adult, 2.5m cable, reusable Pirštinis SpO2 silikoninis minkštas jutiklis, suaugusiems, su 2,5 m kabeliu, daugkartinio naudojimo	1 vnt.	12.2.
01.59.473007	NIBP Tube, reusable (3m) NIBP manžetės prijungimo žarnelė, daugkartinio naudojimo (3m)	1 vnt	12.4.
01.15.040225	2-pin Adult skin temperature probe, reusable (10K/25° C) 2 kontaktų suaugusiųjų odos temperatūros matavimo daviklis, daugkartinio naudojimo (10K/25° C)	1 vnt	12.5.
01.57.471465	ECG cable, 5-lead, clip, Defib, IEC, 3.4m, reusable EKG kabelis, 5 laidų, 5 elektrodams prijungti, Defib, IEC, 3,4 m, daugkartinio naudojimo	1 vnt	12.1.
01.21.064380	Rechargeable Lithium-Ion Battery (10.8V, 2550mAh) (kraunamas Li-jonų akumulatorius (10,8V, 255mAh))	1 vnt.	13.3.