



## needle free valve sterile

att.no.3 pt.2 subpt.2.3

revision 4

date: 2019-03-05

TD-48

amendment of: 2023-02-01

**MANUFACTURER:** ZARYS International Group sp. z o.o. sp.k.

**ADDRESS:** ul. Pod Borem 18, 41-808 Zabrze, Poland, phone +48 32 376 07 00-02, fax. +48 32 370 38 08

### INDICATION:

- intended for closing infusion ports
- device intended for multiple, needle-free injections (liquid drug administration, blood sample taking)
- protects infusion lines from contamination and infection during infusion therapy
- does not require using safety caps
- dedicated for standard luer and luer lock tips



### FEATURES:

- housing made of transparent material - polycarbonate
- high-quality silicone membrane allows multiple port injections
- flat injection surface - smooth and evenly ended port membrane that facilitates disinfection
- metal-free, suitable for mri examination
- male connector of the transparent port is secured with a cap allowing for aseptic removal of the item from the packaging
- latex-free
- DEHP-free
- can be used for 7 days or 350 activations
- shelf life: 3 years
- priming volume: 0,09ml
- high flow rate: 350ml/min
- withstands injection fluid pressure of 3 bars (44 PSI)
- withstands back pressure of 2 bars (29 PSI)
- EO sterilized
- individual packaging: paper-foil

### REFERENCE

| CATALOGUE NUMBER | TYPE        | INTERMEDIATE PACKAGING | BULK PACKAGING |
|------------------|-------------|------------------------|----------------|
| PB-N             | blue        | 100 pcs                | 20 x 100 pcs   |
| PB-T             | transparent | 100 pcs                | 20 x 100 pcs   |