

ENDOPATH® DEXTRUS™

THE MINIMALLY INVASIVE ACCESS SYSTEM
THAT MAXIMISES LAPAROSCOPIC SURGERY.

AN OPEN-AND-SHUT CASE FOR
EXCELLENT PERFORMANCE

SEAL CAP ASSEMBLY:

- Facilitates extracorporeal procedure
- Provides versatility of removing seal cap or specimen through Rapid-adjust Iris Seal

RAPID-ADJUST IRIS SEAL:


- Allows surgeon to custom-adjust iris for comfort and pneumo seal
- Low profile
- Requires fewer device interactions versus the LAP DISC
- Enhanced durability—increased tear strength versus the LAP DISC



ENDOPATH®
dextrus™
MINIMALLY INVASIVE ACCESS SYSTEM

REACH HIGHER

ACHIEVE THE GOALS YOU
THOUGHT WERE OUT OF REACH



AVAILABLE IN 3 SIZES

- Maximise patient selection

STRAIGHTFORWARD SET UP, INSTALLATION AND REMOVAL

- Optimises efficiency

INCREASED TEAR STRENGTH VERSUS THE LAP DISC

- Enhances durability, saves cost

INDICATED FOR USE DURING MINIMALLY INVASIVE AND OPEN SURGERY

- Provides additional surgical options and resources

ENDOPATH® DEXTRUS™ MINIMALLY INVASIVE ACCESS SYSTEM ORDERING INFORMATION

DESCRIPTION	PRODUCT CODE	INSTRUMENTS PER SALES UNIT
Seal Cap Assembly with Accessories	HAP02	3
Fixed-Length Access Retractors		
• Small (up to 4cm)	FLR01	6
• Medium (4-7cm)	FLR02	6
• Large (>7cm)	FLR03	6

©2007. All rights reserved Ethicon Endo-Surgery.