



The new ORTHO VISION™ Analyzer is the first-ever system with Responsive Automation.

## Dynamic Workflow Management

**Efficient:** Designed to respond to your unique workflow, the adaptability to accommodate unpredictable test volume for batch testing.

### CONTINUOUS ACCESS

Load on the fly enables efficient test processing which allows for prioritization of STAT requests with high test volume

TS 7 p.

Ekonominio naudingumo lentelė 6 punktas

### EXTENSIVE MENU

Ability to automate the full range of transfusion medicine testing, including serial dilutions for titration studies and selected cell panels, extended antigen typing and user defined tests

### SELF-REGULATING

Minimum 4-hour walk-away time under fully-loaded conditions and up to 24 hours of productivity allow valuable resources to be dedicated to value-added tasks

### Geresnis naudojimas

Naudojant specialius atidarymo įrankius, analizatoriumi ORTHO VISION galima atidaryti tik tuos kasetės šulinėlius, kurie iš tikrųjų bus naudojami; mažinant atliekų kiekį, išvengiant perkėlimo ir išlaikant aukščiausią saugos lygį.

### IMPROVED UTILITIZATION

Through the use of dedicated punch tools, the ORTHO VISION™ analyzer has the capability to open only those cassette wells that will actually be used; reducing waste while avoiding carry over and maintaining the highest level of safety.

### INTUITIVE INTERFACE

User-friendly design makes it easy to find precisely what you need, which boosts the speed and precision of decision-making.

### ORTHO BIOVUE® CASSETTES

Utilization of trusted column agglutination technology known for its premium accuracy and standardized results offer the required sensitivity and specificity