

# HEALTHPORT® POWER PLP



- Chamber specifically designed to be used for High-flow Injections of contrast media
- Large self-sealing septum facilitates location during skin palpation and reduces risks of puncture concentration
- Transparent and radiopaque elastomeric connector allows rapid, safe and reversible connection of the catheter
- Low profile plastic chamber, also suitable for patients with thick subcutaneous fat layer

## Features

Port

<b>Material:</b>	Titanium + Polysulfone
<b>Shape and dimensions:</b>	circular - base diameter 27 mm
<b>Height:</b>	10.6 mm
<b>Weight:</b>	5.9 g
<b>Internal volume:</b>	0.36 ml

Septum

<b>Material:</b>	Silicone
<b>Diameter available for puncture:</b>	13 mm
<b>Maximum number of punctures:</b>	2000 with 22G Huber needle 1000 with 19G Huber needle
<b>Max Flow Rate:</b>	5ml/s with Max pressure 300PSI, Max Viscosity 12cP

Catheter with rounded tip, marked with a dot every centimeter and a number every 5 centimeters.

Code	Product configuration POWER PLP	Height mm	Weight g	Material	ID mm	OD mm	FR	Length mm	Introducer Dimension FR	Flow mL/min 22G	Flow mL/min 19G
201.24.12.30	All-in-One	10.6	5.9	silicone	1.30	2.70	8	600	9	7	21
201.24.12.20	All-in-One	10.6	5.9	PUR	1.35	2.00	6	600	6.5	8	22
201.24.12.33	All-in-One	10.6	5.9	PUR	1.60	2.70	8	600	9	8	24
201.24.17.30	All-in-One PRO	10.6	5.9	silicone	1.30	2.70	8	600	9	7	21
201.24.17.20	All-in-One PRO	10.6	5.9	PUR	1.35	2.00	6	600	6.5	8	22
201.24.17.33	All-in-One PRO	10.6	5.9	PUR	1.60	2.70	8	600	9	8	24

### All-in-One

#### Package configuration

Port
Catheter
Connector
Straight Huber needle 22G x 1"
Luer connector
Peel-away introducer-dilator
18G introducer needle
0,89 mm J-tip guidewire
Stainless steel tunneler

### All-in-One PRO

#### Package configuration

Port
Catheter
2 connectors
Straight Huber needle 22G x 1"
Luer connector
Peel-away introducer-dilator
18G introducer needle
0,89 mm J-tip guidewire
Stainless steel tunneler
Vein lifter
90° Huber needle with extensions 20G x 1"
10 ml syringe
Blunt tip needle

## Venous Introducer Set

Code		Colour Code
201.51.10.21	Introducer Set 6.5 FR	Red
201.51.10.31	Introducer set 9 FR	Blue

- Peel-away introducer-dilator
- 18G introducer needle
- 0,89 mm J-tip guidewire
- Stainless steel tunneler

žymėjimai ant pakuotės



Reorder Number  
Code article  
Bestellnummer  
Numero di codice  
Número para pedidos  
Número de encomenda  
Numer katalogowy  
Номер заказа



Attention, see instructions for use  
Attention, voir le mode d'emploi  
Achtung! Gebrauchsanweisung beachten  
Attenzione, leggere attentamente le istruzioni per l'uso  
Atención, consulte las instrucciones de uso  
Atenção, consulte as instruções de utilização  
Pozor, prostudujte návod k použití  
Внимание, см. Инструкцию по применению



Use by  
Date limite d'utilisation  
Verfallsdatum  
Data di scadenza  
Utilizar antes de  
Usar até  
Použitelné do  
использовать до

vienkartinis (pažymėta simboliu)



Single use. Do not resterilize  
Produit à usage unique. Ne pas restériliser  
Zum einmaligen Gebrauch. Nicht reesterilisieren  
Monouso. Non resterilizzare  
De un solo uso. No reesterilizar  
Utilização única. Não re-esterilize  
Pro jedno použití. Neprovádějte reesterilizaci  
Одноразового использования.  
Не повторная стерилизация

sterilus (simbolis ant pakuotės)

STERILE EO

Sterilized by ethylene oxide  
Stérilisé à l'oxyde d'éthylène  
Sterilisiert mit Äthylenoxid  
Sterilizzato con gas ossido di etilene  
Esterilizado con óxido de etileno  
Esterilizado por óxido de etileno  
Sterilizováno etylenoxidem  
Стерилизовано оксидом этилена

LOT

Lot number  
Numéro du lot  
Chargenr  
Numero di lotto  
Número de lote  
Número de Lote  
Číslo šarže  
большое число

be latekso (simbolis ant pakuotės)

Not made with  
natural rubber latex

Non fabriqué avec du latex de caoutchouc naturel  
Nicht aus Naturkautschuklatex hergestellt  
Non realizzato con lattice di gomma naturale  
No hecho con látex de caucho natural  
Não fabricado com látex de borracha natural  
Nie je vyrobené z prírodného latexu  
Не сделано с натуральным каучуковым латексом



Manufacturer  
Fabricant  
Hersteller  
Fabbricante  
Fabricante  
Fabricante  
Výrobce  
производитель



MRI Conditional