



Thrombuster™ II

< E version >


6Fr
7Fr

Thrombuster™ II

< E version >

Features

-
- “Big” and “round” suction lumen
 - Minimum flow resistance enhances the suction performance
 - Shortened wire lumen length enhances flexibility of the distal shaft
 - Hydrophilic coating (30 cm from the tip) assures smooth insertion into the blood vessel
 - Removable core wire enhances pushability and anti-kinking feature





Product specifications

Catalogue Number	Guide Catheter Minimum ID Compatibility	Maximum Guide Wire Outer Diameter	Shaft Diameter				Guidewire Lumen Part		Catheter Effective Length
			Outer		Inner		Maximum Outer Diameter	Length	
			Distal	Proximal	Distal	Proximal			
T2R6SA4E	0.070inch (1.78mm)	0.014inch (0.36mm)	3.9Fr (1.30mm)	3.9Fr (1.30mm)	3.0Fr (1.00mm)	3.3Fr (1.10mm)	5.1Fr (1.70mm)	10mm	1400mm
T2R7SA4E	0.080inch (2.03mm)	0.014inch (0.36mm)	4.5Fr (1.50mm)	4.6Fr (1.52mm)	3.6Fr (1.20mm)	3.9Fr (1.32mm)	5.7Fr (1.90mm)	10mm	1400mm

※intended for use in European Territory

- Do not re-use
- Sterilized using ethylene oxide
- Consult instructions for use



Manufacturer

KANEKA CORPORATION

3-2-4, Nakanoshima, Kita-ku, Osaka, Japan

CONTACT TEL. +81-3-5461-3472 FAX. +81-3-5461-3081

EC-Representative

KANEKA PHARMA EUROPE N.V.

Triomflaan 173, B-1160 Brussels, Belgium
TEL. +32-(0)-2-663-0340 FAX. +32-(0)-2-663-0361

