

NIV Interfaces

***FlexiFit™* 431 NIV Full Face Mask** Closed or Open Loop Non-Invasive Ventilation Therapy



14.

The ***FlexiFit™* 431 NIV** Full Face Mask with Under Chin Design builds on the foundation of the unique features common to the ***FlexiFit™*** Series.

The ***FlexiFit™*** Technology has been incorporated into the contoured Silicone Seal to allow a large range of nasal bridge movement and simple one step auto adjustment fitting. The unique ***Glider™*** Strap provides freedom of movement while maintaining a seal. The crown strap and Under Chin Design provide extra stability and maintain a fit even with jaw movement. The Easy-Release Strap allows simple removal of the mask and the “Fit and Forget” Headgear ensures no on-going adjustments are required.

The ***FlexiFit™*431 NIV** Full Face is supplied fully assembled, with all sizes of Silicone Seal in the box.

Product Code	Description	Quantity
HC431 NIV	<i>FlexiFit™</i> 431 NIV Full Face Mask includes: <ul style="list-style-type: none">● Silicone Seal (Small)● Silicone Seal (Medium)● Silicone Seal (Large)● Headgear with Easy Release Cord● Glider™ Strap● Elbows - Non-Vented & Vented● Pressure Port Covers	1 1 1 1 1 2 4

FlexiFit™ 431 NIV Mask Spare Parts

Elbow Kit



Product Code	Description	Quantity
400HC212	Elbows 431 NIV Non-Vented Elbow & Swivel, Vented Elbow & Vented Non-Rebreathing Valve & Instruction Manual (Australia Block)	1 of Each
400HC217	Elbows 431 NIV Non-Vented Elbow & Swivel, Vented Elbow & Instruction Manual (Australia Only)	1 of Each

Silicone Seal




Product Code	Description	Quantity
400HC112	Silicone Seal FlexiFit™ 431 Small	Each
400HC104	Silicone Seal FlexiFit™ 431 Medium	Each
400HC103	Silicone Seal FlexiFit™ 431 Large	Each

Silicone Forehead Pads



Product Code	Description	Quantity
900HC423	Silicone Forehead Pads - fits FlexiFit™ 407 NIV, 431 NIV	2 Pack

Pressure Port Covers



Product Code	Description	Quantity
400HC213	Pressure Port Covers FlexiFit™ 431 NIV	2 Pack

Mask Headgear



Product Code	Description	Quantity
400HC302	FlexiFit™ 431 Headgear	Each

Warning: The components in this catalogue are intended for the Non-Invasive Ventilation mask series. Mixing these parts with any other system may not achieve optimal treatment for a patient using that system. Only use NIV parts for the correct NIV mask.