



The new ORTHO VISION™ Analyzer is the first-ever system with Responsive Automation.

Dynamic Workflow Management

Efficient: Designed to respond to your unique workflow, the adaptability to accommodate unpredictable test volumes for batch testing.

CONTINUOUS ACCESS

Load on the fly enables efficient test processing which allows for prioritization of STAT requests while handling full test volume

TS 7 p.

EXTENSIVE MENU

Ability to automate the full range of transfusion medicine testing, including serial dilutions for titration studies and selected cell panels, extended antigen typing and user defined tests

SELF-REGULATING

Minimum 4-hour walk-away time under fully-loaded conditions and up to 24 hours of productivity allow valuable resources to be dedicated to value-added tasks

Geresnis naudojimas

Naudojant specialius atidarymo įrankius, analizatoriumi ORTHO VISION galima atidaryti tik tuos kasetės šulinėlius, kurie iš tikrųjų bus naudojami; mažinant atliekų kiekį, išvengiant perkėlimo ir išlaikant aukščiausią saugos lygį.

IMPROVED UTILITIZATION

Through the use of dedicated punch tools, the ORTHO VISION™ analyzer has the capability to open only those cassette wells that will actually be used; reducing waste while avoiding carry over and maintaining the highest level of safety.

INTUITIVE INTERFACE

User-friendly design makes it easy to find precisely what you need, which boosts the speed and precision of decision-making.

ORTHO BIOVUE® CASSETTES

Utilization of trusted column agglutination technology known for its premium accuracy and standardized results offer the required sensitivity and specificity