

Hemoclip & reloadable handle Eco Clip

STERILE EO



➤ DISTRIBUTOR

LIFE PARTNERS EUROPE 161 avenue Gallieni – 93170 Bagnole Tel.: +33 1 49 88 01 01 / Fax: +33 1 48 59 70 72 E-mail: life@lifeurope.com Site Internet : www.lifeurope.com	Quality Representative Quality manager Tel: +33 1 49 88 01 01 Fax: +33 1 49 88 83 45 E-mail : quality@lifeurope.com
--	---

➤ PRODUCT INFORMATION

Commercial name:

Eco Clip - Hemoclip & reloadable handle

- Hemostasis clip of 90° or 135° angles jaws, available in a protective cassette. The clip can be easily placed on the reusable handle before being released on the mucosa.
 - True rotary capacity, torque1 for 1.
 - Clip designed with a specific «∞» shape.
 - Comfortable handle with perfect control feeling.
 - In normal use, you can release 100 to 120 clips with one handle*.
- * *Not guaranteed, depending on use and cleaning.*



Indications: this product is intended use for clamping the mucosa tissue of alimentary cavity under the guide of endoscope, in different situations:

- To make marking under the endoscope procedure;
- For hemostasis: - Mucosa/sub-mucosa defects less than 3 cm in diameter; - Bleeding ulcer; - Polyp site is less than 1.5cm in diameter; - Diverticulum in the colon;

References	Maximal opening (mm)	Jaws angle	Jaws length
------------	----------------------	------------	-------------

Hemoclip - Single use

NLS/RHC-8-90	≥ 8	90°	13.5 mm
NLS/RHC-8-135	≥ 8	135°	14.2 mm



Order unit = MDQ (Minimum Delivered Quantity): **30**

Type of individual packaging: **individual blister**

Overall packaging: **cardboard box**

Reference	Catheter length (cm)	Catheter diameter (mm)
-----------	----------------------	------------------------

Reloadable handle - Reusable

NLS/RHCFD-28-230	230	≤ 2.75
------------------	-----	--------



Order unit = MDQ (Minimum Delivered Quantity): **1**

Type of individual packaging: **individual pouch**

Overall packaging: **cardboard box**

➤ **TECHNICAL INFORMATIONS**

Materials:

Handle: ABS
 Clip: Stainless steel
 Clip assembly pack: LCP

Contains:

Latex: Yes No
 Product of animal origin: Yes No
 PVC: Yes No
 Phthalate and DEHP: Yes No
 Nanoparticles: Yes No
 Bisphenol A: Yes No



Sterility:

Sterile: Yes No
 Sterilisation mode: Ethylene oxide gas
 Number of use: Clip: single use



Handle: reusable after sterilization (ultrasonic cleaner, detergent solution, autoclave, Cf. Instructions for use)

Conservation:

Conditions of conservation: The storage environment temperature air, well ventilated , do not store with toxic and dangerous matters
 Use by date: 3 years

Other data:

GMDN class: 34050
 Bar code: GS1-128 HIBC No
 With the information: Reference Lot number Use by date

Use:

Correct and safe use: Cf. Instructions for use
 User advice (instructions for use, precautions and contraindications): Cf. Instructions for use
 Indications (CE mark destination): Cf. Instructions for use

Manufacturer:

Manufacturer: Shangxian Minimal Invasive INC.
 Country: China

CE mark:

CE mark granted by TÜV Rheinland, Certifying body number: 0197
 According to European directive 93/42/CEE – annex II excluding section 4
 Class IIa medical device