

SKINTACT® Monitoring ECG Electrodes for Adults

SKINTACT Monitoring Electrodes offer total versatility at affordable prices. 4 different backings, 2 types of standard gels, a standard and an enhanced adhesive - the SKINTACT range has an electrode to meet all your particular needs in short or long term monitoring.

Protective Foam Backings of short term monitoring electrodes are impermeable to liquids present in surgical procedures.

Breathable Microporous- or Soft Cloth Backings of long term monitoring electrodes provide patient comfort for long term use.

Aqua-Wet Gel ensures immediate signal pick up as it promptly moistures top skin layer.

Aqua-Tac Gel is not only conductive but also adhesive and can be removed without residue.

Easy to open pouch and a handling aid on all SKINTACT monitoring electrodes facilitate use.



FS-521



F-601



F-55



F-TB1



T-60



T-601



W-60



W-601



CT-601

Code	Gel	Backing, Adhesive	Dimensions	Suggested Use
FS-521	Aqua-Tac	Foam, Enhanced	35 x 50 mm	Short Term
F-601	Aqua-Tac	Foam	50 mm	Short Term
F-55	Aqua-Wet	Foam	55 mm	Short Term
F-TB1	Aqua-Tac	Foam	50 x 55 mm	Short Term
T-60	Aqua-Wet	Microporous Tape	60 mm	Long Term
T-601	Aqua-Tac	Microporous Tape	50 mm	Long Term
W-60	Aqua-Wet	Soft Cloth	60 mm	Long Term
W-601	Aqua-Tac	Soft Cloth	50 mm	Long Term
CT-601	Aqua-Tac	Clear Tape	50 mm	Short Term

All SKINTACT® Monitoring ECG electrodes meet the ANSI/AAMI standard EC12 for disposable ECG electrodes.

Packaging Standard:

Aqua-Wet: 30 / pouch • 1200 / box • 4800 / case

Aqua-Tac: 30 / pouch • 1500 / box • 6000 / case