

## DEFIBRINATED SHEEP BLOOD



**Description:** Sheep blood allows detection of hemolytic reactions and supplies the X factor (heme) necessary for the growth of many pathogenic species.

**INTENDED USE:** to enriched culture media

**PRINCIPLES OF THE PROCEDURE:** Sheep blood provides growth factors required for growth of microorganism.

**DIRECTIONS FOR PREPARATION FROM DEHYDRATED PRODUCT:** Add aseptically 50-70 ml of defibrinated sheep blood to sterile medium cooled to 45 - 50°C.

**STORAGE:** bottles – 2-8°C;

**PACKAGING:** cat No. 1000 bottles (50; 100; 200; 500 ml);

**EXPIRATION:** bottles – 21 days;



Graso  
Oddział produkcyjny  
83-211 Jabłowo, ul. Leśna 1 Owidz  
Dział Obsługi Klienta: (58) 562 30 21, fax: (58) 56 279 87  
[zamowienia@graso.com.pl](mailto:zamowienia@graso.com.pl); [www.podloza.pl](http://www.podloza.pl), [www.grasobiotech.pl](http://www.grasobiotech.pl)