


# BIOPSY FORCEPS



## 4. Description

- 3. Fenestred jaws,
- 6. With smooth or notched jaws,
- With or without spike,
- Marking on the catheter for the 160 and 230 cm length
- Biopsy forceps are available with curette from the lot number 2018010101

1.

Minimum delivery quantity (MDQ):  
\* see chart below  
Primary packaging: pouch  
Secondary packaging: Cardboard box



References	Biopsy forceps:	Jaws	Forceps Maxiaml Opening (mm)	Needle	4. Coating	2. Catheter Length (cm)	5. Diameter (mm)	MDQ*
AN/AMHBFA	Tracheo-bronchial	= Standard	5	No	Yes	110	• 1.8	10
AN/AMHBFA-WC	Tracheo-bronchial	= Standard	5	No	No	110	• 1.8	10
1 poz. AN/AMHBFC	Tracheo-bronchial	≡ Alligator	5	No	Yes	110	• 1.8	10
AN/AMHBFC-WC	Tracheo-bronchial	≡ Alligator	5	No	No	110	• 1.8	10
NLS/BFO-18-160	Paediatric	= Standard	5	No	Yes	160	• 1.8	10
NLS/BFOS-18-160**	Paediatric	= Standard	5	Yes	Yes	160	• 1.8	10
NLS/BFO-24-160	Gastric	= Standard	7	No	Yes	160	• 2.4	100
NLS/BFOS-24-160	Gastric	= Standard	7	Yes	Yes	160	• 2.4	100
NLS/BFC-24-160	Gastric	≡ Alligator	5.8	No	Yes	160	• 2.4	100
NLS/BFCS-24-160	Gastric	≡ Alligator	8	Yes	Yes	160	• 2.4	100
NLS/BFO-24-230	Gastric and Colonic	= Standard	7	No	Yes	230	• 2.4	100
NLS/BFOS-24-230	Gastric and Colonic	= Standard	7	Yes	Yes	230	• 2.4	100
NLS/BFC-24-230	Gastric and Colonic	≡ Alligator	5.8	No	Yes	230	• 2.4	100
NLS/BFCS-24-230	Gastric and Colonic	≡ Alligator	8	Yes	Yes	230	• 2.4	100

\*\* Delivery delay on request

➤ **INDICATIONS**

Biopsy forceps are intended for use in collecting tissue samples for histological examination.

➤ **TECHNICAL INFORMATION**

**Materials:**

Jaws: Stainless steel  
 Catheter: Stainless steel  
 Handle: ABS  
 Catheter covering: PE  
 Protecting cap: PVC

**Contains:**

Latex:	<input type="checkbox"/> Yes <input checked="" type="checkbox"/> No	Bisphenol A:	<input type="checkbox"/> Yes <input checked="" type="checkbox"/> No
Product of animal origin:	<input type="checkbox"/> Yes <input checked="" type="checkbox"/> No	Nanoparticles:	<input type="checkbox"/> Yes <input checked="" type="checkbox"/> No
PVC:	<input type="checkbox"/> Yes <input checked="" type="checkbox"/> No	Lead:	<input type="checkbox"/> Yes <input checked="" type="checkbox"/> No
Phthalate and DEHP:	<input type="checkbox"/> Yes <input checked="" type="checkbox"/> No		



**Sterility:**

Sterile: Yes No  
 Sterilization mode: Ethylene oxide gaz  
 Number of use: Single use



**Storage:**

Storage conditions: Store at room temperature  
 Use by date: 2 years

**Other data:**

GMDN code: 38711  
 Barcode: No GS1-128 standard HIBC standard  
 With the informations: Reference Lot number Use by date

**Use:**

Correct and safe use: Cf. Instructions for use  
 User advices: Instructions for use, precautions and contraindications - Cf. Instructions for use



**Manufacturer:**

Manufacturer: ANREI MEDICAL (HZ) Co., LTD  
 Country: People's Republic of China

**CE Mark:**

Notified body: TÜV SÜD, Certifying body N°: 0123  
 Standard: European directive 93/42/EEC  
 Rule: Class IIa medical device, according to rule 6

**IMPORTER**

**LIFE PARTNERS EUROPE**  
 161 avenue Gallieni – 93170 Bagnole, France  
 Phone: + 33 1 49 88 01 01  
 Website: [www.lifeurope.com](http://www.lifeurope.com)

**Order:**

Fax: + 33 1 48 59 70 72  
 E-mail: [life@lifeurope.com](mailto:life@lifeurope.com)

**Vigilance:**

Fax: + 33 1 49 88 83 45  
 E-mail: [quality@lifeurope.com](mailto:quality@lifeurope.com)