




# TKBAG<sup>®</sup>

## SPECIMEN RETRIEVAL BAG



### Safety

Extra-durable bag prevents leaks and tears. Soft ring design reduces tissue damage risks.

### Convenience: Easy 4 U.

Easy in, easy out, easy open, easy close; Saves time!  
Multiple sizes including extra-large

### Pouch with string is optional

Audinių šalinimo sistema, Tinkamas naudoti su 10 mm.

### Quality

Strict, extensive raw material selection and 31 quality inspection steps on every TKBAG.

The TKBAG/ Retriever Bag is used in all types of laparoscopic surgery for removing excised tissues/specimen through a small surgical incision. It can be introduced through 10mm or larger size trocar canula.

It comes in a variety of sizes, (as shown in the table below) and is suitable for a wide range of procedures including, Gallbladders, Appendices, Ectopic pregnancies, Lymph nodes, Ovaries, Lung, Išsiskleidimų o diametras, Maišo ilgis, tons of bowel, Splenectomies, Nephrectomies and more.

Product Code	Bag Diam (mm)	Bag length (mm)	Volume (ml)	Qty/Box
HSD-60	60	120	75	10
HSD-80	80	175	200	10
HSD-80/C	80	175	200	10
HSD-100	100	170	350	10
HSD-130	130	175	700	10



CE 0197  
ISO9001  
ISO13485