

## EXTERNAL DRAINAGE SYSTEM 3

- 82-1730 CODMAN External Drainage System 3 - Complete (sterile) with clear ventricular catheter
- 82-1731 CODMAN External Drainage System 3 - Complete (sterile) without ventricular catheter
- 82-1732 CODMAN External Drainage System 3 Collections bags (5 units)
- 82-1733 CODMAN External Drainage System 3 Leveling Device\*
- 82-1735 CODMAN External Drainage System 3 clear CSF ventricular catheter (1.5mm inner Lumen diameter)
- 82-1738 CODMAN External Drainage System 3 with CODMAN Lumbar Catheter Kit II
- 82-1739 CODMAN External Drainage System 3 clear CSF ventricular catheter (1.9mm inner Lumen diameter)

\*Class 2 laser product



## CODMAN BACTISEAL™ EVD CATHETER

The Anti-microbial Technology Demanded for CSF Drainage. A technically innovative anti-microbial impregnated catheter for external CSF drainage which reduces gram positive bacteria on catheter tube surfaces.

Treated with the patented BACTISEAL Advanced Impregnation Process

- Reducing the potential for bacterial colonization on inner lumen and exterior catheter wall
- Soft, flexible, pliable silicone catheters
- Depth markings for accurate placement
- Requires no technique change
- Available with either normal or large drainage lumen
- Available 'clear' or opaque'

**82-1745** BACTISEAL EVD Catheter Set  
Depth Markings 3cm to 15cm  
Inner Lumen Diameter 1.5mm

**82-1749** BACTISEAL EVD Catheter Set  
Depth Markings 3cm to 15cm  
Inner Lumen Diameter 1.9mm

**82-1750** CLEAR BACTISAL EVD Catheter Set  
Depth Markings 3cm to 15cm  
Inner Lumen Diameter 1.9mm



**82-1745** (Diameter: 1.5mm)



**82-1749** (Diameter: 1.9mm)



**82-1750** (Clear; Diameter: 1.9mm)

### LUMBAR DRAINAGE CATHETER KIT II

The CODMAN Lumbar Drainage Catheter Kit II consists of a lumbar catheter 80cm in length with a closed tip and side holes to allow for drainage without clogging. A Teflon-coated precoiled steel guide wire with an adjustable stop is also included. This kit includes a 14-gauge touhy needle with a clear plastic flash back chamber, centimeter markings, and wings for leverage and ease of use. A female luer lock connector and suture clips made of soft silicone allow for easy connection and attachment to patient.

80cm Radiopaque Catheter

Drainage hole pattern with enlarged hole diameter

100cm guide wire with dispenser

14 gauge touhy needle with 1cm depth markings

Low profile female luer lock and luer lock cap

- 82-1707** CODMAN Lumbar Drainage Catheter Kit II
- 82-1709** CODMAN Lumbar Drainage Catheter Kit II with EDSII (82-1721)
- 82-1738** CODMAN Lumbar Drainage Catheter Kit II with EDSIII (82-1731)

